SÝMBIOZ *Led*

USER MANUAL

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SOPRO COMEG

SOPRO COMEG

USER MANUAL	5
MANUEL D'UTILISATION	29
MANUAL DE USUARIO	53
MANUALE PER UTILIZZATORE	77
BETRIEBSANLEITUNG.....	101
GEbruIKSAANWIJZING	125
INSTRUKCJA OBSŁUGI.....	149

SOPRO COMEG

ENGLISH

SOPRO COMEG

CONTENTS

1	FOREWORD	8
2	ABOUT THE APPLIANCE	9
3	SAFETY INSTRUCTIONS	10
4	REGULATORY INFORMATION	12
5	INSTALLING THE APPLIANCE	13
6	OPERATING PRINCIPLE	15
7	DESCRIPTION OF THE FRONT AND BACK SIDES	16
8	SUGGESTED PROCEDURE FOR DECONTAMINATION	18
9	AFTER-SALES SERVICE AND MAINTENANCE	19
10	TROUBLESHOOTING	20
11	TECHNICAL CHARACTERISTICS	21
12	ELECTRO-MAGNETIC COMPATIBILITY	22
13	SYMBOLS	26

1 FOREWORD

We thank you for your purchase of the SYMBIOZ LED by Sopro-Comeg.

To get the best out of the device while taking all of the necessary precautions, it is essential that you review this manual.

The sentences with with the symbol  are points requiring particular attention.

The sentences beginning with the symbol  are information.

To facilitate the installation and use of the appliance, we wanted to make the manual easier to follow. There are references to the presentation page on the back of the front cover (in the form of [D1] for example) to view the parts of the product concerned more easily.

This user manual is an integral part of the appliance. It must be made available to the user. Correct usage and proper handling of the appliance is made possible by following this manual. The user bears the sole responsibility for damages that result from misuse.

2 ABOUT THE APPLIANCE

SYMBIOZ LED is a LED light source intended to be used by doctors during surgical endoscopy or diagnostic procedures. It is fitted with an electronic iris with fibre optic detection, an automatic thermal protection system, a powerful level memory and an «iris» mode to shut off the light emission.

Its ease of use, its power of lighting and manual fine adjustment of the light intensity make it the ideal multidisciplinary medical tool.

This appliance includes:

- A mains cable
- A Storz type light source adapter
- An Olympus type light source adapter
- A user manual
- A Quickstart guide

Optional accessories:

- Universal light cable
- Adapters for universal light cable, source side (types Storz and Olympus) and endoscope side (Storz, Olympus, Wolf, S&N, ACMI, Stryker types)

We advise you to choose the light cable diameter in the following way:

∅ = 4.8mm for endoscopes of diameter > 4mm

∅ = 3.5mm for endoscopes of diameter ≤ 4mm

 For the United States and Canada, use the supplied "Hospital Grade" cable and connect it to a "Hospital Grade" mains socket.

 This equipment was delivered to you in cardboard packaging. It should be kept for future transport.

 For the United States, Federal law restricts the sale of this product to physicians or under their order.

3 SAFETY INSTRUCTIONS

Read the user manual.

- The appliance should only be connected to a supply network with a grounding system.
- Observe the conditions of use and storage.
- The appliance must be opened exclusively by a competent technician authorized by the manufacturer.
- Do not insert any metal objects into the appliance to avoid any risk of electrocution, fire, short circuit or dangerous emissions.
- Do not expose the appliance to splashes of water or store in a damp place.
- Only use the accessories delivered with the appliance or offered as an option by the manufacturer.
- Do not place any heavy objects on the appliance.
- This unit is not sterile.
- If the power cord is damaged, immediately turn the appliance off. It is dangerous to operate this appliance with a damaged cable.
- To disconnect the cable, pull it out by the plug. Never pull on the cable itself.
- Disconnect the appliance from the mains if you do not plan to use for a few days or more.
- In the area of installation, the presence of high interference fields (strong magnetic fields for example) should be avoided. These interference fields can result in malfunctions of the light source.
- Be sure to provide sufficient air circulation to avoid overheating inside the appliance: at least 15 cm around the appliance. Do not cover the side and back air vents and ensure that the legs of the device are there.
- Use only the disinfectant methods recommended in section 8.
- Before each use, make sure that the appliance has no rough surface, cutting edge, or protuberances which could cause problems for safety.
- The surface temperature of the light cable ends can exceed 41°C (after a few minutes of use). Therefore, avoid holding area in contact with the patient.
- To avoid any risk of electric shock, this appliance must be connected only to a grounded mains supply socket.
- This appliance may not be modified without authorization of the manufacturer. - If the medical device is modified, an inspection and a test must be performed to ensure that the medical device complies with the safety instructions.
- Check the compatibility of your endoscope with the manufacturer before use.
- This appliance must be used on individuals (patients) suitable for undergoing a procedure of endoscopy;
- Never look at the light output or the end of the light cable.
- Do not insert anything other than a light cable into the housing [C1]. Doing so may damage the optical system.
- The use of other cables or accessories than those specified may cause an increase in emissions or decreased immunity of the appliance.
- After use of the appliance, when the fibre is removed from its light-guide, the temperature on the metal coupler of the fibre is very high and can cause burns.
- The light power at the output of the fibre can cause eye injuries. Caution when handling the

3 SAFETY INSTRUCTIONS

fibre when the appliance is in use.

- Do not place the distal end of the light cable or of the endoscope directly on the patient or on any other inflammable matter (linen, gauze, operative field, etc.) because it can be very hot and cause burns.
- A multisocket or extension cable must not be connected to the EM system.
- **The light source is intended for professional use in an operating theatre.**
- **Before you connect the appliance, check whether the mains voltage and frequency indicated on the appliance correspond with the values of the mains power supply.**
- **This light source can be connected only to the Sopro camera.**
- **Do not drop the appliance. If the appliance is dropped, do not reconnect it but return it to your authorized distributor or directly to the after-sales department.**

 **After use, be sure to remove the light cable carefully without touching the distal tip which can be very hot and cause burns.**

 It is advisable to have a second light source available in the operating room to use if an absence or degradation of performance is observed.

 This product is equipped with Group 1 LEDs according to the standard IEC 62471. Do not look directly at the light to avoid any risk of eye damage.

The potential equalization plug located behind the medical device may be used to equalize the earth potential of the medical device with that of all the devices connected to the mains in the environment. Use the common earth of the hospital or the building.

Connect the potential equalisation cable to the equipotential C3 plug on the back of the appliance.

 The devices that connect to the inputs/outputs must comply with IEC 60950-1

4 REGULATORY INFORMATION

4.1 COMPLIANCE

This appliance was designed and made by a company with a certified quality system. It meets the requirements of European directive 93/42/EEC concerning medical devices. Therefore, it particularly meets electrical safety standards (IEC) and Electromagnetic Compatibility (EMC) ad hoc.

4.2 ELECTROMAGNETIC INTERFERENCE AND ELECTROSTATIC DISCHARGE

Although this product complies with EMC standards, it is possible that it may disrupts other devices in very specific circumstances, or itself be disturbed by other devices or an unfavourable electromagnetic environment.

To avoid these situations, the following is recommended:

- Ensure the quality of mains electricity (and particularly the grounding of all appliances and trolleys).
- Keep the appliance away from electromagnetic sources (eg: compressor, motor, transformer, HF generator, etc.).

4.3 MATERIOVIGILANCE

Like any medical device, this device is subject to the provisions of materiovigilance. Any serious malfunction must therefore be reported to the competent authorities and the manufacturer in the shortest possible time and with the highest possible accuracy. Manufacturer contact Information: refer to the last page of the manual.

4.4 END OF LIFE

This unit bears the recycling symbol in accordance with directive 2002/96/EEC concerning Waste Electrical and Electronic Equipment (WEEE or WEEE).

By disposing of this device properly you will help to prevent any harmful consequences to the environment and to human health. The symbol  on the appliance or the accompanying documents indicates that this product may never be treated as household refuse. It must, consequently, be taken to a refuse collection centre charged with recycling electrical and electronic equipment. For disposal, observe the standards relating to the dumping of refuse in force in the country of installation.

For more details about dealing with, recovering and recycling this product, please contact your nearest dealer who will tell you what to do.

5 INSTALLATION OF THE APPLIANCE

The purpose of this device is the lighting of the work area through an endoscope.

- This medical device is intended to be used by a qualified surgeon during endoscopic explorations. Its installation requires no special training. Please refer to the instructions in this manual;
- Install the unit on a stable surface and take the various accessories essential for its operation;
- If you use a unit, ensure that it leaves at least 15 cm of air around the appliance.
- Connect the power cord to the socket of the appliance **[C2]** ;
- Put switch **[S2]**, on the 1 position. The device is in standby mode **(the green light associated with the Standby button [S1]) should be flashing, and a logo appears on the touch screen (L1))** ;
- Press the Standby button **[S1]** **or on the touch screen** to start the appliance. The green light comes on and a logo indicating that a light cable should be inserted appears on the screen **[L1]**. Pressing the Standby button again will put the appliance in the standby position.

Remark: If the light source is connected to the Sopro camera through the communication cable supplied on the socket (C4), starting the latter will automatically start the light source. The latter may however be put back into standby mode, or switched on independently.

Remark: If the light source is connected to the Sopro camera, the triggering of the white balance from the Sopro camera will start the light source automatically. The latter may however be put back into standby mode, or switched on independently.

- Push the light cable into the socket for this purpose **[C1]**, you should hear a locking «click». **The logo asking for a light cable to be inserted will disappear to leave room for the main screen of the interface. The LED is switched off at this stage. Pressing the «On» button of the interface will cause the LED to deliver a level of intensity equal to the last one used.**

Remark: Starting the Sopro camera using a personal SCard will cause the LED to start at the level of intensity recorded by the user (after the «On» button being pressed by the user).

- Connect the other end of the light cable to your endoscope.
Before installing an endoscope (or a fibroscope) in front of the camera lens, check the condition of this endoscope (clean distal lens or eyepiece).
- **Press the «On» button of the interface to exit the iris mode and turn on the light. This button changes state at this stage to become an «Off» button. Press this button to switch off the LED.**

5 INSTALLATION OF THE APPLIANCE

 The light source is originally equipped with a «Storz» type light cable. **An «Olympus» type adapter is available behind the product in its housing [A2].**

 If no light cable is connected to the device when it is running, a logo explaining this information will be displayed on the screen and the LED will be off. Only then connect a light cable to use the appliance.

 **After use, be sure to remove the light cable carefully without touching the distal tip which can be very hot and cause burns.**

6 OPERATING PRINCIPLE

6.0 GENERAL OPERATING PRINCIPLE

An LED is supplied with current according to the intensity setting. This LED emits light which is focused toward the exit of the optical fibre. The optical fibre transmits light towards the endoscope to light the working area.

6.1 ADJUSTING THE LIGHT INTENSITY

After inserting the light cable, a logo indicating the «iris» mode (LED off) will be displayed on the screen [L1].

Pressing the «On» button «Turn On» or clicking on the logo makes it possible to leave this mode and deliver the light. The logo, initially empty in «iris» mode will fill accordingly, indicating the adjusted level of light intensity.

If the Sopro camera is connected to this light source through the supplied communication cable, the triggering of the white balance from the camera will switch on the LED, if it had been turned off. The movement of the intensity cursor of the interface to the right allows the increase of the light intensity. To the left, it decreases the intensity.

The empty logo indicates that the source gives a minimum intensity while the complete logo indicates that it gives a maximum light intensity.

 The LED takes about 10 seconds to stabilize its colour temperature. You are advised to wait for this period before doing the white balance on the camera.

 It is also advisable to do a white balance if changing the intensity adjustment.

6.2 AUTOMATIC MODE

When the «automatic control» mode is activated, the light intensity is automatically adjusted according to the lighting conditions of the environment. It remains possible, in this mode, to activate the manual intensity adjustment slider, which will deactivate the automatic mode. In addition, this mode can also be deactivated by switching «OFF» by the corresponding switch.

6.3 POWER LEVEL MEMORY

To put it the LED source in position «iris», simply press the «Off» button or click on the logo. This button will then become an «On» button. Press it to return to the previously set light power. Similarly, when the LED source has been disconnected from its power supply, it will restart on the last level of light intensity used.

7 DESCRIPTION OF THE FRONT AND BACK

7.1 DESCRIPTION OF THE AIR VENTS

The **air vents [A1] and [A3]** on the back and side of the appliance must not be obstructed, to avoid any risk of overheating.

The light source is equipped with an automatic safety device which stops the lighting if the internal temperature becomes excessive. A pictogram showing the activation of this safety device will be displayed on the screen **[L1]**.

7.2 MAINS SOCKET

The electrical power supply of the appliance is via a supplied standard IEC power cord.

7.3 MAINS SOCKET - FUSES

The appliance receives its electrical power supply via the mains socket located on the back **[C2]**, which must be connected to the mains by the supplied power cord. This socket has a fuse slot and the general power switch **[S2]**. When replacing a fuse, it is imperative to disconnect the product from the network and use a fuse of the same type. The T of «T2.5A» means «timed». Use only fuses T2.5AL - 250V - UL.

7.4 INDICATIONS AND SYMBOLS ON THE BACK

The signs and symbols on the back of the product allow you to identify it in accordance international standards IEC 60601-1, IEC 60601-2-18, IEC 60417 and EN 980 (see paragraph 13).

7.5 COMMUNICATION PORT

A communication port **[C4]** is available for the connection of this light source to the Sopro camera.

A service port **[C5]** is for updating the internal software by an authorised Sopro technician.

A1	Back air vent
A2	Connection tip of for Olympus type light cable
A3	Side air vent
C1	Connection tip of for Storz type light cable
C2	Mains socket
C3	Equipotential socket
C4	Connection to the Sopro camera

7

DESCRIPTION OF THE FRONT AND BACK

C5	Department
F	Fuse compartment
L1	Touch-screen
S1	Standby button
S2	Power button
I	Label

8 SUGGESTED PROCEDURE FOR DECONTAMINATION

 The LED source is a medical device, not heat resistant and it cannot withstand immersion. Therefore it is disinfected using a non-woven cloth soaked in a detergent-disinfectant.

 As decontamination is linked to products, methods and/or selected tools, it remains under the full responsibility of the staff concerned.

 This unit cannot be autoclaved.

 The existing alkaline solutions for the disinfection of certain medical devices are PROHIBITED for the disinfection of this device.

 **To ensure hygiene between patients and non-contamination, the light cable must be changed from one patient to another following the rules of endoscopic practice (accessory not supplied which must always be sterilized before each use).**

9 AFTER-SALES SERVICE AND MAINTENANCE

No particular maintenance operation is needed for this appliance.

 **Reminder:** misuse is not covered by the warranty. If a fault persists and it is necessary to return the appliance to the after-sales service, take care to ship it in its original packaging. Similarly, it is appropriate to send the complete device back (controls, cables). Please send a short explanatory note about the fault with the dispatch note.

 The equipment must be disinfected before returning for repair. When returning the equipment, check its condition and make notes on the delivery note if necessary, confirming them with the carrier by registered letter within 48 hours. If equipment shipped by us suffers damage during transport, the amount of the repair is payable by the carrier if the complaints have been made within the deadline, otherwise by the recipient. In the event of an incident, please contact our customer service or our nearest sales representative to you.

10 TROUBLESHOOTING

10.1 THE LED OF THE STANDBY SWITCH [S1] DOES NOT LIGHT UP WHEN SWITCHED ON

- Check that the mains cable is properly connected to the mains and the appliance and the power switch [S2] on the back is in the operating position («I»).
- Check that the fuses are in good condition (use only fuses T2.5A - 250V - UL).

10.2 THE SOURCE LIGHTS UP BUT THE LIGHT FLOW IS INSUFFICIENT

- Check the condition of your light cable and your optics.
- Check that the intensity setting is not too low.
- Verify that the light cable is pressed in fully (it needs to stop clearly when it is fully home).

10.3 THERE IS NO MORE LIGHT BUT THE FANS CONTINUE TO OPERATE

The light source is equipped with a safety device which shuts off the supply to the LED if the temperature in the appliance is too high. A pictogram showing this safety device will be displayed on the screen. Once the temperature has dropped again, the source returns to Standby position and can be used again.

- **Check that the light cable is correctly inserted in its housing.**
- Check that there is enough space around the source to ensure sufficient cooling (15cm all around).
- Check that nothing is blocking the ventilation grilles on the back and below the source.
- The light source is in «iris» mode. Press the **«On» button** or move the intensity adjustment slider to release the light.

If a fault persists and it is necessary to return the appliance to the after-sales service, take care to ship it in its original packaging after having disinfected it.

10.4 THE LIGHT APPEARS/DISAPPEARS INTERMITTENTLY

- Check that the light cable is pressed fully home and that the connector [C1] is properly screwed in.

In the event of a fault, please contact your authorized dealer or the after-sales department for any assistance.

11 TECHNICAL CHARACTERISTICS

Main functions

- LED technology
- Colour temperature: > 6500°K
- LED typical life duration 40000H
- Compatible light cable Storz and Olympus
- Electronic Iris
- Light intensity adjustment: Manual and automatic
- Anti-blinding system with detector of light cable
- Standby mode
- Automatic thermal protection system
- Power level memory
- Supply voltage: 100-230VAC / 50 - 60 Hz
- Electrical consumption: 160 VA
- Dimensions (W x H x D): 310 x 136 x 385 mm
- Weight: 5000g.

Environment

- Operating temperature: +10°C / +40°C
- Operating humidity: 30% to 75 %
- Transport and storage temperature: -10°C / +40°C
- Transport and storage humidity: 20 to 85%
- Operating, storage and transport atmospheric pressure: 800hPa to 1 060hPa
- Not protected against water spillage (IPX0)
- Not suitable for use in the presence of inflammable anaesthetic gas mixture with air, oxygen or nitrous oxide.

Regulations

Class 1, type CF electrical safety to European Directive 93/42/EEC - Class I
Complies with international standards IEC 60601-1 ; IEC 60601-2 -18; IEC 60417 and EN 980. This product is equipped with group 1 LEDs to standard IEC 62471. Do not look directly at the light to avoid any risk of eye damage.

12 ELECTRO-MAGNETIC COMPATIBILITY

12.1. Manufacturer's guide and declaration - electromagnetic emissions		
The equipment referenced as «LED Source» is intended to be used in the electromagnetic environment specified below. The user must ensure that it is used in this environment.		
Emissions test	Conformity	Electromagnetic environment - guide
RF emissions CISPR 11	Group 1	This «LED Source» product is used in the RF energy only for its sub-systems. It therefore emits very low RF energy and it is not probable that it interferes with nearby electronic equipment. This «LED Source» product is suitable for use in all places other than domestic and those directly connected to the public low-voltage power supply network supplying buildings for domestic use.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complying	

12 ELECTRO-MAGNETIC COMPATIBILITY

12.2. Manufacturer's guide and declaration - electromagnetic Immunity			
The equipment referenced "LED Source" has been designed to be used in the electromagnetic environment specified below. The user must ensure that it is actually used in this environment.			
Immunity test	CEI 60601 Severity level	Compliance tests	Electromagnetic environment - guide
Electrostatic discharge IEC 61000-4-2	±6 kV via Contact ±8 kV via the air	± 6 kV ± 8 kV	The floor should be wood, concrete or ceramic tiles. If the ground is covered with synthetic material, the relative humidity should be at least 30%.
Fast transient peaks IEC 61000-4-4	±2 kV electrical lines ±1 kV input/output	± 2 kV ± 1 kV	The quality of the main power supply should be that of a typical commercial or hospital environment.
Electric shock IEC 61000-4-5	Differential mode ±1 kV Common mode ±2 kV	± 1 kV ± 2 kV	The quality of the main power supply should be that of a typical commercial or hospital environment.
Electrical problems, short cuts of current and variations of voltage IEC 61000-4-11	<ul style="list-style-type: none"> • <5% UT – for 10 ms • 40% UT – for 100 ms • 70% UT – for 500 ms • <5% UT – for 5 s 	<ul style="list-style-type: none"> <5% UT 10 ms <40% UT 100 ms <70% UT 500 ms <5% UT 5 s 	The quality of the electrical power supply must be that of a typical commercial or hospital environment. If the user of this product must be able to continue to work during power cuts, it is recommended that this product be powered by UPS or battery.
Magnetic Field on the network frequency (50/60 Hz) IEC61000-4-8	3 A/m	3 A/m	The magnetic field at the frequency of the network frequency of the network must be at a level characteristic of a place (50 / 60 Hz) in a typical commercial or hospital environment.
Remark: UT is the nominal value of the electrical voltage applied during the test.			

12 ELECTRO-MAGNETIC COMPATIBILITY

Manufacturer's guide and declaration - electromagnetic safety			
The equipment referenced "LED Source" is intended to be used in the electromagnetic environment specified below. The user must ensure that it is used in this environment.			
Safety test	IEC 60601 Severity level	Compliance tests	Electromagnetic environment - guide
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	<p>The portable and mobile RF communications equipment should not be used at a distance of this product - including cables - lower than what is recommended, calculated by applying the formula that corresponds with the frequency of the transmitter.</p> $d = 1,16 \sqrt{P}$ $d = 1,16 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2,33 \sqrt{P} \text{ 800 MHz to 2.5 GHz}$ <p>In which P is the maximum power output from the transmitter in watts (W), assigned by its manufacturer and (d) is the recommended separation distance in meters (m).</p> <p>The levels of fields emitted from transmitters RF fixed - which must be established by an electromagnetic measurement in situ - must be less than the compliance level in each frequency band. Interference may occur with devices on which is placed.</p> <p>The following symbol: </p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	
<p>Remark 1: At 80 MHz and 800 MHz, it is appropriate to use the highest frequency band.</p> <p>Remark 2: These recommendations may not be applicable to all situations. The propagation of electromagnetic waves is modified by the absorption and reflection due to the structures, objects, and people.</p> <p>For transmitters whose maximum output power is not listed on the table above, the recommended separation distance d, in metres (m) can be established using the equation applicable to the frequency of the transmitter, in which P is the maximum output power of the transmitter in watts (W) assigned by the manufacturer of the transmitter.</p>			

12 ELECTRO-MAGNETIC COMPATIBILITY

12.3 Recommended distances between reportable and mobile communication systems and this product

The equipment referenced "LED Source" has been designed to be used in an electromagnetic environment in which the RF interference emitted is controlled. The user of this device can help to prevent electromagnetic interference by keeping a minimum distance between the of portable and mobile RF communications systems (transmitters) and this device, as recommended below, as a function of the maximum output power of the communications system.

Maximum output power assigned of the transmitter in W	Separation distance as a function the transmitter frequency		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1,16 \sqrt{P}$	$d = 1,16 \sqrt{P}$	$d = 2,33 \sqrt{P}$
0.01	0.116	0.116	0.233
0.1	0.366	0.366	0.736
1	1.16	1.16	2.33
10	3.66	3.66	7.36
100	11.6	11.6	23.3


Remark 1: At 80 MHz and 800 MHz, it is appropriate to use the highest frequency band.

Remark 2: These recommendations may not be applicable to all situations. The propagation of electromagnetic waves is modified by the absorption and reflection due to the structures, objects, and people. For transmitters whose maximum output power is not listed on the table above, the recommended separation distance d , in metres (m) can be established using the equation applicable to the frequency of the transmitter, in which P is the maximum output power of the transmitter in watts (W) assigned by the manufacturer of the transmitter.

13 SYMBOLS

	S	Button
	L	Indicator light
	C	Socket
	I	Label
		Manufacture date
		Manufacturer
Class 1		Class I product
		Conforming to European directive 93/42/EEC
		An equipotential earth socket conductor other than a protective earth or a neutral conductor, allowing direct connection between the electrical equipment and the equalizing bar of the installation potential. Please consult standard IEC 60601-1 3rd edition.
	T	Timed Fuses UL/CSA Device of CF type
		CF-type device
		Electronic and electrical equipment put on the market after 13/08/2005. This symbol indicates that this product must not be processed with household waste.
		Read the user manual.

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