SÝMBIOZ *Cam*

USER MANUAL
MANUEL D'UTILISATION
MANUAL DE USUARIO
MANUALE PER UTILIZZATORE
BETRIEBSANLEITUNG
GEBRUIKSAANWIJZING
INSTRUKCJA OBSŁUGI

SOPRO COMEG

SOPRO COMEG

USER MANUAL	5
MANUEL D'UTILISATION	37
MANUAL DE USUARIO	69
MANUALE PER UTILIZZATORE	101
BETRIEBSANLEITUNG.....	133
GEBRUIKSAANWIJZING	165
INSTRUKCJA OBSŁUGI.....	197

ENGLISH

SOPRO COMEG

CONTENTS

1	PREFACE	8
2	INTRODUCTION TO THE APPARATUS	9
3	SAFETY INSTRUCTIONS	10
4	REGULATORY ADVICE	12
5	INSTALLING THE APPARATUS	13
6	OPERATING GUIDELINES	16
7	DESCRIPTION OF THE FRONT AND BACK PANELS	22
8	SUGGESTED DECONTAMINATION PROCEDURES	24
9	AFTER-SALES SERVICE AND MAINTENANCE	26
10	TROUBLESHOOTING	27
11	TECHNICAL CHARACTERISTICS	29
12	ELECTROMAGNETIC COMPATIBILITY	30
13	SYMBOLS	34

1 PREFACE

We thank you for your purchase of the SYMBIOZ CAM by Sopro-Comeg.

In order to make the best use of it while having all the necessary precautions at your disposal, it is essential for you to become acquainted with this manual.

The sentences displaying the  symbol correspond to points requiring special attention.

The sentences displaying the  symbol represent items of information.

To facilitate installation and use of the apparatus, we have attempted to make the device manual more practical. Consequently, references to the product presentation page on the inside cover (like [D1] for example) will be provided for easier viewing of the relevant product parts.

This user manual is an integral part of the device. It must be made available to the user. Proper use and correct handling of the device entail following these instructions. The user alone shall be responsible for any damage that may result from improper use

2 INTRODUCTION TO THE DEVICE

This is a 3Cmos FULL HD medical camera.

Its ergonomically designed sensor, automatic shutter, good sensitivity, excellent resolution and faithful colour-rendering make it ideal as a medical tool for rigid or flexible endoscopy.

This camera comes with:

- A control unit
- A camera head with integrated 15-34 mm zoom lens
- Video cables
- A power cable
- A peripheral cable
- Two communication (SOPROCOM) cables
- A set of 6 Smartcards
- A user manual
- A Quickstart guide

 **For the United States and Canada, the “hospital grade” power cord supplied and the connector must be used with a “hospital grade” mains socket.**

 This equipment has been delivered to you in packaging to be retained for use if transporting the device.

3 SAFETY INSTRUCTIONS

 Read the User Manual.

 This user manual is an integral part of the apparatus. It must be made available to the user. Proper use and correct handling of the apparatus entail following these instructions. The user alone shall be responsible for any damage that may result from improper use.

- Comply with the conditions of use and storage;
- Do not expose the apparatus to excessive dust;
- The apparatus should be opened only by a competent technician authorized by the manufacturer;
- Do not insert metal objects into the apparatus. This is to avoid any risk of electrical shock, fire, short-circuit or hazardous emissions;
- Do not expose the apparatus to splashed water or to damp;
- Use only the accessories supplied with the apparatus or recommended as options by the manufacturer;
- This camera is not a flame-proof device. Therefore, it must not be used in the presence of flammable anaesthetics;
- This apparatus was not designed to be used in an ionising environment;
- This apparatus is not sterile;
- The surface temperature of the camera can reach 41°C (after a few minutes' use). Therefore, avoid any contact with this area.
- The camera is intended for use by qualified, trained professionals, in an operating theatre.
- Prior to plugging in the apparatus, check that the mains voltage and mains frequency indicated on the apparatus correspond to the power system values.
- Prior to activation, ensure that neither the apparatus nor the power cable is damaged. Deteriorated cables and connectors must be replaced immediately.
- Do not drop the apparatus. If the apparatus falls, do not reconnect the apparatus but send it back to your authorised distributor or directly to the Sopro after-sales service department.
- No additional multiple-socket outlets or extension cords must be connected to the EM system.

 To avoid any risk of electrical shock, this device must be connected only to a power system equipped with protective grounding.

 Apparatuses connecting to the input/output ports must comply with the IEC 60950-1 standard.

 Any modification of this device without authorisation from the manufacturer is prohibited. If the medical device is modified, an inspection and a test must be carried out to ensure that the medical device complies with the safety regulations.

3 SAFETY INSTRUCTIONS

 Check with the manufacturer regarding the compatibility of your endoscope and your light source, prior to use.

 It is advisable to have a second surgical camera in the operating theatre so that action can be taken if the device fails to perform or if a deterioration in performance is noticed.

 This device is to be used on individuals (patients) fit to undergo an endoscopic procedure.

The potential equalization plug located at the back of the medical apparatus can be used to equalize the grounding potential of the medical device with that of all the apparatuses plugged into the power main in the environment. Use the shared grounding system in the hospital or the building.

 Apparatuses connecting to the input/output ports must comply with the IEC 60950-1 standard.

4 REGULATORY ADVICE

4.1 COMPLIANCE

This product was designed and manufactured by a company with a certified quality system. It meets the requirements of European directive 93/42/CEE, on medical devices. Consequently, it particularly meets the standards of electrical safety (IEC) and electromagnetic compatibility (CEM) ad hoc.

4.2 ELECTROMAGNETIC INTERFERENCES AND ELECTROSTATIC DISCHARGE

Although this product complies with CEM standards, it may in very special circumstances interfere with other devices, or itself be the object of interference from other apparatuses or an unfavourable electromagnetic environment.

Image transfers from the camera to a monitor may be susceptible to electromagnetic disturbances and electric system disturbances.

In order to avoid these situations, it is advisable to:

- Ensure of the quality of the electric power system (especially the grounding of all apparatuses and medical carts).
- Keep the apparatus away from electromagnetic sources (e.g. compressors, motors, transformers, HF generators, etc.).

4.3 MEDICAL DEVICE VIGILANCE

Like any medical device, this apparatus is subject to the stipulations governing medical device vigilance, and therefore any serious malfunction must be reported to the competent authorities and to the manufacturer as quickly and as accurately as possible. For manufacturer contact details, refer to the last page of the manual.

4.4 END OF LIFE CYCLE

This apparatus carries the recycling symbol in compliance with European directive 2002/96/CEE on Waste Electrical and Electronic Equipment (DEEE or WEEE). By correctly disposing of this apparatus you are helping to prevent harmful effects on the environment and on human health.

The  symbol displayed on the apparatus and on the accompanying documentation indicates that this product cannot under any circumstances be treated as household waste. It must therefore be delivered to a waste collection centre for the recycling of electrical and electronic equipment.

In disposing of it, please comply with the waste elimination norms in effect in the country where it is installed.

For further details on the treatment, recovery and recycling of this apparatus, kindly contact your nearest retailer who will advise you on the procedure to follow.

5 INSTALLING THE APPARATUS

This medical device is intended for use by a qualified surgeon during endoscopic exploration. No special training is required to install it. Kindly refer to the instructions in this manual.

Place the control unit on its base and have to hand the various accessories required to operate it.

5.1 CONNECTIONS

- Connect the power cord to the mains socket of the camera [C2] and a video lead to one of the video output ports: We recommend selecting the DVI output [C6] for connecting the main surgical monitor.
- Connect the second end of the video lead to the video input port on the monitor.
- If you wish to control certain functionalities of the Sopro insufflator or the Sopro light source, connect these to the SOPROCOM socket on the camera's control unit [C9] intended for this purpose, using the two cables provided.

Connecting the endoscope to the sensor

- Turn the ring on the endoscope-holder on the camera head clockwise.
- Once the endoscope is inserted, release the ring.
- If the sensor has to be locked into position on the endoscope, the ring should be turned in the opposite direction.
- Connect the potential equalization cable to the equipotential plug C13 at the back of the apparatus.

5.2 ACTIVATION

 To optimize the quality of the images obtained with this product, always be sure to connect the camera head before turning on the control unit.

- Connect the power cable to the plug [C2] at the back of the camera.
- Connect the DVI-D cable to the corresponding output port [C6] at the back of the product and to the corresponding input port on the monitor.
 - ⇒ The HDS DI output produces a progressive 1080 signal
- If needed, connect other monitors or peripherals, preferably using the HDS DI output ports [C8]
 - ⇒ HDS DI output ports produce an interfaced 1080 signal
 - ⇒ HDS DI output ports are preferred when a very long cable is necessary.

5 INSTALLING THE APPARATUS

- Turn the activation switch [C2] to position "1"
 - ⇒ The camera takes about 30 seconds to come on. The Sopro Comeg logo is displayed on the screen during this interval.
 - ⇒ The flashing LED standby button, and the appearance of "standby" on the screen indicate that the camera is properly activated
- Connect the camera head to the control unit via the dedicated connector [C1]
- Touch the screen, OR press the Standby button to access the main menu.
- Select a surgical speciality from Laparoscopy, Cystoscopy, Hysteroscopy, ENT, Arthroscopy and Flexible Endoscopy, by clicking the associated symbol.
 - ⇒ Once the selection is made, the camera is in surgery mode
- After connecting the light cable and the endoscope, film a white surface and launch the White Balance by pressing the dedicated button on the touch screen interface, or the camera head's pre-programmed button
 - ⇒ When the surgical monitor appears indicating "AWB OK" and the chronometer starts, this signals that the white balance has successfully been set.
 - ⇒ The chronometer shows the operating time by default (the starting point being when the white balance is set). However, this chronometer can be set beforehand if necessary, independently of the white balance.
- The video settings have been defined for each surgical speciality. However, they can be modified, if necessary, during surgery.
 - ⇒ 3 of these (Brightness, Sharpness and Red Gain) can be directly accessed from the main screen. All the others (Gamma, Blue Gain, Colour, Red Phase, Gain, Noise Reduction, and Aperture) can be accessed by clicking the "More Settings" button. This is a drop-down list. To scroll through the parameters, touch the left section (the parameter headings) with your finger.
- Press "Close settings" to exit and return to the main screen.
- You can at any time during surgery to return to the factory settings for the surgical speciality by clicking the "Specialities presets" button.
- Press the "STOP" button at the end of the surgery.

You can at this stage:

Return to the main menu by clicking the "Completed" button;

OR

Create a user profile from the settings used during the operation;

OR

Log in if you wish to add this set of settings to your profile.

5 INSTALLING THE APPARATUS

 Before connecting an endoscope (or fibrosopes) to the camera lens, check the condition of this endoscope (clean distal lens, setup if necessary)

5.3 ACTIVATING THE DEVICE FROM A SMARTCARD

- Scan a Smartcard on the read area provided on the lexan housing (next to the Sopro Comeg logo)
- Scanning a standard Smartcard:
 - ⇒ Activate the camera directly on the surgical screen, in the surgical application corresponding to the Smartcard scanned.
- Scanning a personal Smartcard:
 - ⇒ Activate the camera directly on the surgical screen, in the surgical preset corresponding to the Smartcard scanned. If the user profile associated with the Smartcard scanned shows more than one preset, then an intermediate stage for selecting the desired preset will be necessary.
- Scanning an administrator Smartcard:
 - ⇒ Activate the camera in the administrator menu

6 OPERATING GUIDELINES

6.1 SELECTING THE SURGICAL SPECIALITY

This camera can be used in various mini-invasive surgical operations.

Select the speciality from the pre-entered list. This loads the video parameters for the application accordingly:

LAPAROSCOPY
 HYSTEROSCOPY
 CYSTOSCOPY
 ENT
 ARTHROSCOPY
 FLEXIBLY ENDOSCOPY

Note 1: Each of these specialities can be accessed by scanning the Smartcard on its dedicated area [Z1] on the control unit. This can be done at any time. If the Smartcard is scanned while a surgical operation is in progress, then confirmation will be requested.

Note 2: If you have entered customized parameters into a user profile beforehand, then you can select this user profile either by using the "Log in" button, or by scanning your personal Smartcard on the dedicated area [Z1] of the control unit if the control unit is linked to your personal profile.

6.2 CONFIGURING THE FUNCTIONS OF THE CAMERA HEAD BUTTONS

This operation is performed on the corresponding page of the menu interface, which can be accessed from the dedicated "Camera head" button:

Click the button you wish to configure (B1, B2 or B3).

Select the type of press (short press or long press) and the corresponding function desired, and confirm.

Repeat this procedure for each of the three buttons.

The configuration selected for each of the three buttons will be entered into the user's personal profile.

Below is the comprehensive list of the configurable functions:

Type of press	
Short press	Long press
No action	No action
Brightness +	Brightness +
Brightness -	Brightness -
Light source - Intensity +	Light source - Intensity +
Light source - Intensity -	Light source - Intensity -
Peripheral 1	Peripheral 1

6 OPERATING GUIDELINES

Peripheral 2	Peripheral 2
Mirror	Mirror
Freeze frame	Freeze frame
Recording images on the camera (if equipped with an integrated recorder)	AWB
Record / Stop video recording on the camera (if equipped with an integrated recorder)	Insufflators (RUN/High Flow Rate/STOP)
	Light source (ON/OFF)
(a) Recording images on the PC	Recording images on the camera (if equipped with an integrated recorder)
(b) Record / Stop video recording on the PC	Record / Stop video recording on the camera (if equipped with an integrated recorder)
	(a) Recording images on the PC
	(b) Record / Stop video recording on the PC

By default, the following configuration is applied (factory setting):

B1 (short press / long press): Brightness - / Light source

B2 (short press / long press): Brightness - / Insufflator

B3 (short press / long press): No action / AWB

Mirror mode: An "M" is displayed on the monitor to indicate that the mirror mode (inverted image) is active, and disappears when it is deactivated.

6.3 WHITE BALANCE

The white balance can be set in 2 ways.

- From the touch screen interface, click the white balance button.
- From one of the camera head buttons, depending on your programming (only possible with long press)

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- Once the camera is paired with the endoscope, the ambient light reduced to a minimum and the light source active, film a white surface.
- Set the white balance; the BALANCING IN PROGRESS message is displayed on the screen.
- Continue to film the white surface until the BALANCE OK message is displayed on the screen.

Note 1: Setting the white balance causes the chronometer to be activated, signalling the start of the surgery and therefore the availability of additional functionalities, particularly those linked to the image quality settings.

Note 2: If the Sopro light source is connected to the Sopro camera by the communication cable supplied, then setting the white balance will cause the LED to light up if it was previously off.

6 OPERATING GUIDELINES

⚠ In the event that the white balance is set during an operation, i.e. inside the patient, the endoscope must be pulled out of the trocar and the white balance must be set again using gauze.

6.4 FOCUSING

Once the endoscope is connected and the light source is activated, slowly turn the lens focus ring to find a position where the objects being observed are sharp.

i Focusing on a fairly distant point helps obtain adequate field depth for the operation and therefore avoids having to do retouching too regularly.

6.5 FUNCTIONALITIES

6.5.1 VIDEO SETTINGS

Once the white balance is set, or the "Start" button has been pressed, access to the video settings is authorised.

In this way, it is possible to work directly on brightness, sharpness and red gain. Additional settings are available by clicking the "More settings" button.

The following parameters can be modified:

BRIGHTNESS	This parameter varies the general brightness of the image
BLUE GAIN	Colorimetry shift from white to blue
RED GAIN	Colorimetry shift from white to red
COLOUR	Modifies colour saturation
NOISE REDUCTION	Reduces the noise in an image
APERTURE	Analysis aperture on the shutter (electronic iris): image surface on which the camera analyses the quantity of light received, in order to adapt the shutter speed to obtain the proper exposure time. This function is used to obtain good exposure even when the image has areas of widely varying illumination (e.g. a dark orifice with very light surrounding areas requires a small aperture).
GAMMA	Modifies the dynamics of the dark areas: +1 reduces the difference between the light portions and the dark portions, therefore making these dark areas appear less so. -1 does the opposite: the dark portions are made even darker.
SHARPNESS	Modifies the sharpness of the image to obtain a 3D effect or on the other hand, softens the image with the anti-moiré filter for flexible endoscopes.
GAIN	Modifies the brightness level in cases of mediocre lighting conditions
RED PHASE	

Click the “Close settings” button once to return to the previous menu.

6.5.2 OTHER PARAMETERS

HD Recording

If your camera is equipped with an integrated HD recorder, you can set HD image / video recording from the interface menu using the corresponding buttons or from the buttons on the camera head if these have been programmed accordingly.

Insufflators

If your camera is connected to the Symbioz Flow insufflator via the communication cable provided, you can select from the parameters menu whether or not to display information about the operating status of the insufflators on the surgical monitor’s OSD.

In particular:

Feedback on instantaneous pressure in the pneumoperitoneum;
Feedback on warning messages from the insufflators.

In addition, one of the buttons on the camera head can be configured (long press only) to control the Symbioz Flow insufflators according to the following sequence: Start insufflation in low flow / Change to high flow / Stop insufflation.

6.6 COMPLETION

Once the operation is completed, press the “STOP” button.

You will be given the option:

To create a user profile;

To log in, if you have not already done so;

To save any modifications to the settings in your personal profile, if you have not already done so;

To exit.

Note:

If you are already logged in, you can opt to save the modifications to your personal profile.

In particular you will be able to save:

- Modifications to the video settings;
- Modifications to the Sopro light source intensity setting (last setting used);
- Modifications to the Sopro insufflator pressure setpoint (last setting used, with 15 mm Hg limitation, and mode);

6 OPERATING GUIDELINES

- Modifications to the settings on the surgical monitor supported models only);
Once the selection has been made, you can overwrite the existing preset, or create another one.

6.7 GENERAL SETTINGS MENU

The general settings menu is used to:

- Set the light intensity on the touch screen
- Select the display messages and the display of the instantaneous pressure in the pneumoperitoneum, from the Sopro insufflator.
- Select the language

6.8 INFORMATION MENU

The information menu is used to access a rapid start-up product guide.

6.9 CONNECTIONS MENU

The connections menu is used to obtain the connection status of all of the system's devices. In particular, the camera head, the Sopro light source, the Sopro insufflators, the Sopro PC Panel and the Surgical Monitor.

Note: The icon for this menu changes to orange when at least one of the products on the system is disconnected.

6.10 CREATING A USER PROFILE

- Click the "Log in or create a profile" button
- Click the "Create a profile" button
⇒ Activating the camera with a blank Smartcard gives direct access to the following stage
- Click the text field and enter the desired profile name
- Click the "Return" button
- Select a surgical speciality. This will form the basis for the video parameters.
- Click the "Create a profile" button OR scan a blank Smartcard on the area provided for this purpose on the lexan housing (next to the Sopro Comeg logo)
- A message is displayed for a few seconds indicating that the device supports the creation of a profile.

6 OPERATING GUIDELINES

- You can then:

- ⇒ Rename the profile and/or the preset by clicking the "Manage" button
- ⇒ Delete one (of the) profile preset(s) by clicking "Manage"
- ⇒ Use the preset created by clicking on it
- ⇒ Log out by clicking the "Log out" button

Notes:

If a blank Smartcard has been scanned at least once (during activation and/or while saving a created profile), then the profile (and the presets) created will be automatically linked to this Smartcard. A padlock icon will be displayed on this profile. This means that:

- ⇒ The presets locked in this profile cannot be modified without the corresponding Smartcard.
- ⇒ The names of the presets and the profile name cannot be modified without the corresponding Smartcard.

7 DESCRIPTION OF THE FRONT AND BACK PANELS

7.1 MAINS SOCKET

The camera receives electrical power by means of the mains plug located on the back panel [C2], which must be connected to the mains power supply by the lead supplied with the camera. This plug carries a fuse trap as well as a master switch for turning on the apparatus. The camera uses T2.5AL - 250V fuses.

When replacing a fuse, the camera must be unplugged and the same type of fuse must be used. The T in "T2.5A" means "time-delay". Only use fuses marked UL/CSA.

7.2 VIDEO OUTPUT PORTS

The following video output ports are available on this camera: DVI [C6], HD-SDI [C8], RGBS [C5], Y/C [C10] and composite [C11].

i The DVI output port provides the best image quality in 1080p format over short distances. **i** HD-SDI output ports are used to transport the video signal over long distances (screen at a distance of up to 50 m...)

7.3 PERIPHERAL OUTPUT PORTS

Two (2) peripheral output ports are available [C12] and these are used with a 3.5 mm jack lead for managing two peripherals (for example: a printer and a video recorder) from the buttons on the camera head.

⚠ The apparatuses connecting to the "VIDEO OUT" and "PERIPHERAL" ports must be compliant with standard IEC 60950.

7.4 KEYS AND SYMBOLS ON THE BACK PANEL

The keys and symbols on the back panel are for identifying the camera in accordance with international standards IEC 60601-1, IEC 60601-2-18, IEC 60417 and EN 980 (See Appendix 1)

The table below describes the back and front components

C1	Connector for the Camera head
C2	Main power socket
C3	RS232 output port (connection to Panasonic EJ-MLA26E surgical monitor)
C4	Service line
C5	RGBS video output port
C6	DVI-D video output port
C7	Equipotential socket
C8	Ethernet output port (Connection to the Sopro insufflators)

7 DESCRIPTION OF THE FRONT AND BACK PANELS

C9	Connections to Sopro products (light source and insufflators)
C10	Y / C video output port
C11	Composite video output port
C12	Peripheral connections
C13	Equipotential plug
L1	Touch screen
S1	Standby button
S2	Power button
Z1	Read area for Smartcard

8 SUGGESTED DECONTAMINATION PROCEDURES

Disinfecting the camera control unit

 The camera control unit is a non-heat-resistant medical device and does not withstand immersion; consequently, it should be disinfected using a non-woven medium saturated with a disinfectant detergent.

 This apparatus is not drip-proof. Do not splash liquid on this product.

 As decontamination is linked to the products, methods and/or tools used, decontamination is the sole responsibility of the staff concerned.

 The apparatus must always be disconnected before any cleaning procedure is carried out.

 This apparatus is not autoclavable.

 Existing alkaline solutions for the disinfection of certain medical devices are NOT RECOMMENDED for disinfecting this apparatus.

Disinfecting the camera head

Pre-disinfection cleaning	Disinfection		
Immersion in an enzymatic cleaning solution (Aniosyme DD1, Hexanios G + R or surface cleaning ((Anios wipes))).	Immersion in a glutaraldehyde solution (Anios Laboratory - Steranios 2%)	Sterilisation using the STERRAD process	Sterilisation with ethylene oxide

 Any other disinfection methods are prohibited, and the manufacturer accepts no liability for any damage caused by using such methods.

 It is advisable to use a disposable, sterile, protective cover on the camera regardless of the sterilisation procedure followed and throughout the surgical operation.

 The camera is not autoclavable.

 The camera is not compatible with automatic washer-disinfectors.

8 SUGGESTED DECONTAMINATION PROCEDURES

 Existing alkaline solutions for the pre-disinfection of certain medical devices are NOT RECOMMENDED for pre-disinfecting our cameras.

 It is essential for the parts coming into contact with the disinfectant to be thoroughly rinsed.

 Use soft non-woven cloths for wiping the lenses dry so as not to scratch them.

 The procedures described in this section are given as guidance, and cannot under any circumstances be substituted for official recommendations or directives.

 The connector on the camera equipped with its sealing cap can be immersed. Traces of water are visible just as the connector seal is removed (runoff and condensation). It is therefore essential to dry the inside of the cap well before reusing it.

Disinfecting the lens surfaces

A dirty lens surface can interfere with observation. The lens surface must not have any marks or smudges. To avoid scratching the lens surface, never use abrasive cloths or sponges to clean it.

To remove dust, dirt and other impurities not originating from the patient, wipe clean with a soft, lint-free cloth moistened with ethanol or 70% isopropyl alcohol. When cleaning, wear gloves for protection from chemical products.

Ensure that the equipment is completely dry before use.

 Use soft non-woven cloths to clean the lenses so as to avoid scratching them.

9 AFTER-SALES SERVICE AND MAINTENANCE

This device does not require any special maintenance operations.

Note: misuse of the device is not covered by the warranty.

If a fault persists and the device has to be returned to the after-sales service department, ensure that it is shipped in its original packaging.

Similarly, it is advisable to return the device in its entirety (control unit, power cables).

Kindly attach to the shipping order a short explanatory note about the fault noted.

The equipment must be disinfected before it is returned for repair. When returning the equipment, check its condition and note any exceptions on the shipping order if necessary, and confirm these with the shipper by registered mail within 48 hours. If any equipment shipped by us suffers damage during transport, the total amount for repairs will be charged either to the shipper if notification of the exceptions has been given within the deadline, or failing such, to the recipient. In the event of an incident, kindly contact our customer service department or your nearest sales representative.

10 TROUBLESHOOTING

10.1 THE INDICATOR LIGHT ON THE STANDBY SWITCH [S1] DOES NOT LIGHT UP ON START-UP

- Check that the mains socket at the back of the camera is properly connected to the power supply and that the general switch at the back of the camera is in the operating position ("I").
- Check that the fuses are in good condition (only use the fuses specified on the back panel).

10.2 THE INDICATOR LIGHT ON THE STANDBY BUTTON [S1] COMES ON BUT NO IMAGE APPEARS ON THE SCREEN

- Check that the sensor is properly connected to the control unit (if not, a colour bar will be displayed), and that it is filming a workable image.
- Check that the control unit on the camera is correctly connected to the monitor (cable in good condition, plugs correctly pushed in).
- If the video channel has several components, connect the camera directly to the monitor to verify that the problem is not with one of the peripherals.
- Check that the monitor is in fact switched on, wired to the correct video input and that the image settings are at the mid-point position (colour, light and contrast).
- Check the light source, the light cable and the endoscope.

10.3 THE IMAGE IS BLURRED AND WHITE ALL OVER

- Check that there is no condensation on the lens.
- Check that the sensor is not aimed at an object that is too bright.
- Check the lens setup.

i If there is no image or if the image is not displayed, unplug and re-plug your apparatus or reset the camera to the factory settings.

10.4 THE IMAGE IS EXTREMELY LIGHT OR DARK

- Check that the BRIGHTNESS parameter is not at the maximum or minimum setting.
- Check that the APERTURE used is not too wide or too narrow.
- Check that the camera head is correctly connected to the control unit.
- Check that the light intensity from the light source is not too low or too high.
- Turn off and turn back on the control unit.

If the fault persists and it becomes necessary to return the camera to the after-sales service department, be sure to ship it in its original packaging after disinfecting it. Similarly, it is advisable to return the camera in its entirety (control unit, sensor, lens, cables). Kindly attach to the shipping order an explanation concerning the fault noted.

10 TROUBLESHOOTING

 The equipment must be disinfected before it is returned for repair.

When returning the equipment, check its condition and note any exceptions on the shipping order if necessary, and confirm these with the shipper by registered mail within 48 hours. If any equipment shipped by us suffers damage during transport, the total amount for repairs will be charged to the shipper if notification of the exceptions has been given within the deadline, or otherwise to the recipient.

11 TECHNICAL CHARACTERISTICS

CONTROL UNIT

Main functions

- Class 1 device, CF-type applied part
- 3 programmable buttons
- Application and user MENUS
- Automatic compensation function for defective pixels
- RF transmitter:
 - Type : RFID
 - Frequency : 13,56 MHz
 - Modulation feature : OOK
 - Transmitter power : 70mW (+18 dBm)

Video output

- CVBS: 1.0 V [P-P] / 75 ohms Negative Sync,
- Y/C: Y: 1.0 V [P-P] / 75 ohms, C: 0.286 V [P-P] (NTSC), C:300 Vp-p V [P-P] (PAL) / 75 ohms Negative Sync,
- RGB: 3x0.7 V [P-P] & 0.3V [P-P] / 75 ohms
- HD-SDI: 0.8 V [P-P] / 75 ohms SMPTE292M
- DVI-D 3x0.7Vcc/& 0,.3Vcc/75Ω

Environment

- Operating temperature: +10°C / +40°C
- Operating humidity: 30% to 75 %
- Transport and storage temperature: -10°C / +45°C
- Transport and storage humidity: 20 to 85 %
- Operating, storage and transport atmospheric pressure: 700hPa to 1 060hPa
- Control unit is not drip-proof (IPXO)
- Not suitable for use in the presence of a flammable anaesthetic mixture with air, oxygen or nitrous oxide
- Compliant with European directive 93/42/CEE.
- Compliant with international standards IEC 60601-1; IEC 601-2-18; IEC 60417 and EN 980
- Dimensions (W x H x D): 310 x 136 x 385 mm
- Control unit weight: 4300g
- Electric power supply: 100 - 230 VAC @ 50 Hz-60 Hz
- Electricity consumption: 38-50 VA
- Two T2.5AL - 250V fuses (Marked UL/CSA only)
- Continuous operation

Camera head

- CF-type applied part
- Sensor: 3 CMOS
- Acquisition in native 16:9 format
- Progressive scan technology
- Resolution: HDTV 1080p format: 1920x1080 pixels
- Definition: 1080 lines
- Sensitivity: 5 lux
- Signal/Noise ratio: 54dB
- Cable length: 3m
- Sensor dimensions (D x H): 149 x 53 mm
- Sensor weight: 270g

12 ELECTROMAGNETIC COMPATIBILITY

Manufacturer's guide and declaration – electromagnetic emissions		
This camera is intended for use in the electromagnetic environment specified below. The user must ensure that it is used in this environment.		
Emission test	Compliance	Electromagnetic environment - Guide
RF emissions CISPR 11	Group 1	This camera only uses radio energy for its internal functions. Consequently, its RF emissions are very weak and are not likely to interfere with nearby electronic equipment.
RF emissions CISPR 11	Class A	This camera must be used on premises other than domestic premises and premises directly connected to a low-voltage public energy distribution network providing power to buildings for domestic purposes.
Harmonic Emissions EN 61000-3-2	Class A	
Voltage fluctuations / Flicker EN 61000-3-3	Compliant	

12 ELECTROMAGNETIC COMPATIBILITY

Manufacturer's guide and declaration – electromagnetic immunity			
This camera is intended for use in the electromagnetic environment specified below. The user must ensure that it is used in this environment.			
Immunity test	IEC 60601 Level of gravity	Compliance level	Electromagnetic environment - guide
Electrostatic discharges EN 61000-4-2	± 6 kV in Contact ± 8 kV in air	± 6 kV ± 8 kV	The floor must be wooden, concrete or tiled. If the floor is covered with a synthetic material, the relative humidity must be at least 30%.
Rapid transient peaks EN 61000-4-4	± 2 kV for power lines ± 1 kV for input/ output lines	± 2 kV ± 1 kV	The quality of the main power supply must be that of a commercial or typical hospital environment.
Electric shocks EN 61000-4-5	Differential mode ± 1 kV Shared mode ± 2 kV	± 1 kV ± 2 kV	The quality of the electric power supply must be that of a commercial or typical hospital environment.
Power dips, power failures, short power interruptions and variations in power voltage EN 61000-4-11	<ul style="list-style-type: none"> • <5% UT – for 10 ms • 40% UT – for 100 ms • 70% UT – for 500 ms • <5% UT – for 5 s 	<ul style="list-style-type: none"> <5% UT 10 ms <40% UT 100 ms <70% UT 500 ms <5% UT 5 s 	The quality of the main power supply must be that of a commercial or typical hospital environment. If the user of this camera requires it to continue functioning during interruptions of the main power supply, it is advisable for this camera to be powered by an inverter or a battery.
System frequency magnetic field (50/60 Hz) IEC61000-4-8	3 A/m	3 A/m	The system frequency magnetic field must be at a level that is characteristic of location in a commercial or typical hospital environment.

12 ELECTROMAGNETIC COMPATIBILITY

RF conducted EN 61000-4-6	3 Vrms 150 kHz at 80 MHz	3V	<p>Portable and mobile RF communication devices must not be used nearer to the equipment, reference "Mono CDD camera", including the cables, than the recommended separation distance, calculated using the applicable formulas based on the transmitter frequency. Recommended separation distance.</p> <p>$d = 1,16 \sqrt{P}$</p> <p>$d = 1,16 \sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 2,33 \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power from the transmitter in Watts (W), assigned by the manufacturer of the transmitter and d is the recommended separation distance in metres (m).</p> <p>The field strengths emitted by fixed RF transmitters, determined from electromagnetic measurements of site a, must be lower in terms of compliance in each frequency band.</p> <p>Interference may be caused by nearby devices displaying the following symbol: </p>
RF emitted EN 61000-4-3	3 V/m 80 MHz at 2.5 GHz	3V/m	

NOTE : Ur is the nominal value of the power voltage applied during the test.

NOTE 1 : at 80 MHz and 800 MHz, the higher frequency band applies.

NOTE 2 : these recommendations may not be applicable in all situations. The propagation of electromagnetic waves is modified by absorption and reflection due to structures, objects and people.

a The field strength of fixed transmitters, such as radio telephone (mobile / wireless) and land mobile radio base stations, amateur radio systems, AM/FM radio communication systems and TV systems cannot in theory be evaluated with accuracy. To obtain an evaluation of the electromagnetic environment due to fixed RF transmitters, measurements must be taken of the site. If a field strength measured in the environment in which this camera is used exceeds the applicable compliance levels shown above, this camera must be checked for proper functioning. If it observed to be functioning abnormally, additional measures must be taken, such as a change of direction or repositioning of the referencing equipment.

b Outside of the 150 kHz to 80 MHz frequency band, the field strength must be lower than 3 V/m.

12 ELECTROMAGNETIC COMPATIBILITY

Recommended separation distances between portable and mobile RF communication systems and this camera			
This camera is intended for use in an electromagnetic environment in which the RF interference produced is controlled. The user of this camera can help avoid electromagnetic interference by maintaining a minimum distance between the portable and mobile RF communication systems (transmitters) and this camera, as recommended below, based on the maximum output strength of the communication equipment.			
Assigned maximum output strength of the transmitter in W	Separation distance depending on the frequency of the transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1,16 \sqrt{P}$	$d = 1,16 \sqrt{P}$	$d = 2,33 \sqrt{P}$
0.01	0.116	0.116	0.233
0.1	0.366	0.366	0.736
1	1.16	1.16	2.33
10	3.66	3.66	7.36
100	11.6	11.6	23.3

For transmitters for which the maximum output strength is not shown in the table above, the recommended separation distance d , in metres (m) can be determined by using the equation applicable to the frequency of the transmitter, where P is the maximum output strength of the transmitter in Watts (W) assigned by the manufacturer of the transmitter.


NOTE 1 : at 80 MHz and 800 MHz, the separation distance given in the higher frequency band applies.

NOTE 2 : these recommendations may not be applicable in all situations. The propagation of electromagnetic waves is modified by absorption and reflection due to structures, objects and people.

13 SYMBOLS

	S	Button
	L	Screen
	C	Socket
	I	Label (under the product)
		Date of manufacture
		Manufacturer
		Protective grounding
Class 1		Class I product
		Compliant with European directive 93/42/CEE
		An equipotential grounding socket conductor other than protective grounding or a neutral conductor, for direct connection between the electrical equipment and the potential equalisation busbar. Kindly refer to standard CEI 60601-1, 3rd edition.
	T	UL/CSA time-delay fuses
		Video output ports and peripheral control output port
		Data input and peripheral output
		CF-type device
		Electronic and electrical equipment marketed after 13/08/2005. This symbol indicates that this product is not to be treated with household waste.
		Read the user manual
		RF transmitter

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