


User manual
Manuel d'utilisation
Manual de usuario
Manuale per utilizzatore
Betriebsanleitung

SOPRO *184*

INTUITIVE
CONNECTION

medical technologies

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● User manual	5
● Manuel d'utilisation	29
● Manual de usuario	53
● Manuale per utilizzatore	77
● Betriebsanleitung	101

English

Contents

1	Foreword	8
2	Presentation of the video camera	9
3	Safety instructions	10
4	Notice required by law	11
5	Installation of the video camera	12
6	Connecting the endoscope	13
7	Operation in automatic mode	14
8	Operation in manual mode	16
9	Description of the back panel	18
10	Procedure for disinfecting the sensor	19
11	After-sales service	20
12	Troubleshooting	21
13	Technical characteristics	22
14	Electromagnetic compatibility	23
15	Symbols	27

1 Foreword

Thank you for the trust you have shown in us by acquiring this device.
To derive maximum benefit from it whilst also taking all the requisite precautions, you must read this manual.

Sentences featuring the  symbol are points that require special attention.

Sentences featuring the  symbol provide information.

To facilitate the installation and use of the video camera, we have done our best to make the manual for it more practical. Consequently, cross-references to the presentation page on the back of the cover ([D1] etc.) are offered so you can view the parts of the product concerned more easily.

2 Presentation of the video camera

The device in question is a fully DIGITAL colour CCD mono video with remote electronics. Its ergonomic compact sensor, its automatic shutter, and its sensitivity and excellent resolution, as well as its outstanding colour rendering, make it an ideal medical tool for rigid or flexible endoscopy.

This video camera includes:

- A control channel unit,
- A mains cable,
- A BNC video cable (composite),
- A Y/C video cable,
- An RGB video cable,
- A USB2 cable,
- A BNC video cable (if there is an SDI option),
- A stereo jack cable (for controlling peripherals),
- A user manual,
- 1 CD

Optional items:

- A 1/4" CCD straight sensor with a C mount and featuring programmable buttons, along with a 3m cable and a connector for the control channel unit and its airtight plug.
- A lens with a focal length of 18 mm, a C mount, with a ring for holding an endoscope.
- 2 SDI outlets.

 **For the United States and Canada, use the Hospital Grade cord and the connector that must be fitted with a Hospital Grade mains plug.**

 **Federal law restricts this device to sale by or on the order of a doctor.**

 This equipment will be delivered to you in a carton or a suitcase, which should be kept should you need to transport this equipment.

3 Safety instructions

Read the user's manual

- Follow the storage and user conditions;
- Do not expose the device to excessively dusty conditions;
- The device should only be opened by a skilled technician authorised by the manufacturer;
- To avoid any risks of electric shock, fire, short-circuiting or hazardous emissions, do not insert any metal objects into the device;
- Do not splash water onto the device or leave it in a location that is too humid;
- Only use the accessories provided with the device or offered on an optional basis by the manufacturer;
- This video camera is not a flammable device. Consequently, it must not be used in the vicinity of inflammable anaesthetics.
- This device was not designed for use in an ionising environment.
- This device is not sterile.
- The video camera's surface temperature can reach 41°C (after a few minutes of use). Consequently, avoid any contact with this area.

 **To avoid any risk of electric shock, this device should only be connected to a power supply network with earthing protection.**

 **Devices connected to the inputs/outputs must comply with the IEC 60950-1 standard.**

 **It is prohibited to modify this device without the manufacturer's authorisation. If the medical device is modified, an inspection and a test must be carried out in order to make sure that the medical device complies with the safety instructions.**

 **Check the compatibility of your endoscope with the manufacturer, along with your light source, and do so prior to any use.**

 It is recommended that you have a second surgical video camera available in the operating unit so that work can be done should there be a failure, or should a deterioration in performance levels be observed.

 **This device must only be used on individuals (patients) fit for undergoing an endoscopy procedure.**

4 Notice required by law

4.1 Compliance

This product was designed and manufactured by a company with a certified quality system. It fulfils the requirements of European Directive 93/42/EEC, relating to medical devices. Consequently, in particular it complies with electrical safety standards (IEC) and ad hoc ElectroMagnetic Compatibility (EMC) standards.

4.2 Electromagnetic interferences and electrostatic discharge

Although this product does comply with EMC standards, it is possible under very special circumstances that it may disrupt other devices, or that it may itself be disrupted by other devices or an unfavourable electromagnetic environment.

In order to avoid these situations, the following are recommended:

- Ensuring the quality of the electrical network (and particularly the earthing of all the devices and carriages)
- Keeping the device away from electromagnetic sources (for example the compressor, motor, transformer, HF generator, etc.).

4.3 Medical device vigilance

As is the case with any medical device, this device is subject to the provisions relating to medical device vigilance; any serious malfunction must therefore be reported to the competent authorities and to the manufacturer as soon as possible, and as accurately as possible.

For the manufacturer's contact details, refer to the last page of the manual.

4.4 End of lifecycle

This device bears the recycling symbol, in accordance with European Directive 2002/96/EEC concerning Waste Electrical and Electronic Equipment (WEEE).

By disposing properly of this device, you will contribute to preventing any harmful consequences for the environment and human health.

The  symbol on the device or on the documentation accompanying it indicates that this product may not be treated like household waste under any circumstances. Consequently, it must be disposed of at a waste collection center responsible for recycling electrical and electronic equipment waste.

For disposal purposes, comply with the standards relating to disposing of waste in force in the country where the device is installed. To obtain more details about the topic of treatment, recovery and recycling of this device, please contact the competent office in your borough, the waste collection company, or contact your retailer directly.

5 Installation of the video camera

This medical device is intended for use by a qualified surgeon as part of endoscopic explorations. Installing it does not require any special training. Please refer to the instructions in this manual. Place the communication channel on a stable support and take its various accessories.

- Connect the power supply cord to the video camera's mains socket [C2].
- Connect a video cable to one of the video outputs: Preferably, use the Y/C output [C8], the SDI output [C3 & C4], or the RGB output [C7], which provide a better result than the composite output [C9].
- Connect the second end of the video cable to the monitor's video input.

Connect the video camera's USB2 output [C5] to the computer's USB2 input.

- Place the switch [S9] in position 1. This will place the video camera in stand-by mode (the green indicator [L1] should then flash).
- Hook up the sensor connector to the front panel of the control channel unit [C3].
- Power up the monitor.
- Press the button [S1] to start up the video camera (the green indicator will light up continuously). Pressing it once more will return it to stand-by position.
- An image will appear on the monitor and a driver installation message will appear on the computer; when it does, select the driver delivered with the video camera.

It is prohibited to carry out a surgical operation with the USB2 output; a video output must be connected to the monitor.

Only the EMPIA2820 driver should be installed for the video camera's USB2 output. (The driver is available with the CD delivered.)

Then it will be possible to check that the video camera is operating properly and in particular to make sure that the monitor's settings are completely satisfactory: when you disconnect the sensor, a bar test pattern will appear (except on the USB outlet [C5]).

Before fitting an endoscope (or a fibrescope) to the front of the camera lens, you should check the condition of the endoscope in question (check that the distal lens or eye-piece is clean).

The functions can be accessed using a single adjustment wheel [S2]. Turning the wheel allows you to browse through the Menu and pressing on the wheel enables you to validate an item.

6 Connecting the endoscope

● 6.1 With the straight sensor

Bring the pegs [S6] on the ring [R1] together and then insert the endoscope.

Once the endoscope is inserted, release the two pegs.

If the sensor's position has to be locked in relation to the endoscope, you should push the pegs apart until there is a clicking sound.

Then all you have to do is bring the pegs together in order to unlock the position of the video camera.

7 Operation in automatic mode

 This video camera is equipped with an automatic mode for adjusting the video parameters. In just two steps, the video camera will be operational. This mode is activated by default.

7.1 Selection of the type of endoscope

Select the rigid endoscope (endoscope with diameters of 1.9 to 10 mm fitted with lenses) or the flexible endoscope (endoscope fitted with fibre optics) in the Optics tab on the screen [L2] by pressing on the wheel [S2].

7.2 Focusing

Use the wheel [S4] on the sensor to carry out the focusing.

Once the endoscope is connected and the source is activated, turn the wheel slowly in order to find a position where the objects observed are in sharp focus.

 Carrying out focusing quite far away enables you to obtain sufficient field depth for the operation, which does away with the need to regularly carry out focusing adjustments.

7.3 White balance

There are two ways of initiating the white balance.

1. From the control channel unit, by pressing the wheel when the AWB white balance icon is selected.
2. From the button [S7] from the straight sensor.

You should carry out the following steps:

- Once the video camera is connected to the endoscope, the ambient light is reduced to a minimum, and the light source is active, film a white surface.
- Activate the white balance; the message "Balance In Progress" will appear on the OLED screen of the control channel unit, or 4 white squares will appear on-screen.
- Continue to film the white surface until the "Balance In Progress" message on the OLED screen of the control channel unit disappears, or until the 4 white squares disappear.
- At the end of the white balance process, a square matching the size of the endoscope used will appear on-screen; this defines the exposure window calculated by the video camera.

 If you use one of the buttons on the head of the video camera, keep the button pressed until the "Balance In Progress" message is obtained and then release it; this time delay will avoid unintentional pressing on the button.

7 Operation in automatic mode

7.4 Tools menu

 Depending on whether you have selected Auto or Manual mode, the contents of the Tool menu will differ.

Select the Tool tab on the screen [L2] in order to gain access to the following Menu:

- Languages: EN: English; FR: French; DE: German; ES: Spanish; IT: Italian.
- Up/Down key: Programming the functions that can be activated from the head of the video camera. In automatic mode, the bottom button [S7] is always programmed for carrying out the white balance [AWB]. It is left to the user to choose one of the following options for parameterisation of the top button: endoscope, peripherals, mirror, no action.

 An "M" will be displayed on the monitor in order to indicate that the Mirror mode has been activated and it will disappear when this mode is deactivated.

- Blue shift: Shift in the white colorimetry towards the blue end of the scale.
- Red shift: Shift in the white colorimetry towards the red end of the scale.

8 Operation in manual mode

8.1 Selecting the manual mode

Select the Mode tab and choose "Manual".

From the list, select the type of endoscope used:

- Rigid 1.9 mm
- Rigid 2.7 mm
- Rigid 4 mm
- Rigid 5 mm
- Rigid 10 mm
- Semi-Rigid
- Flexible
- User's choice

8.2 White balancing

There are two ways of initiating the white balance.

1. From the control channel unit, by pressing the wheel when the AWB white balance icon is selected.
2. From the button [S7] from the straight sensor.

You should carry out the following steps:

- Once the video camera is connected to the endoscope, the ambient light is reduced to a minimum, and the light source is active, film a white surface.
- Activate the white balance; the message "Balance In Progress" will appear on the OLED screen of the control channel unit, or 4 white squares will appear on-screen.
- Continue to film the white surface until the "Balance In Progress" message on the OLED screen of the control channel unit disappears, or until the 4 white squares disappear.
- At the end of the white balance process, a square matching the size of the endoscope used will appear on-screen; this defines the exposure window calculated by the video camera.

 If you use one of the buttons on the head of the video camera, keep the button pressed until the "Balance In Progress" message is obtained and then release it; this time delay will avoid unintentional pressing on the button.

8 Operation in manual mode

8.3 Tools menu

 Depending on whether you have selected Auto or Manual mode, the contents of the Tool menu will differ.

Select the Tool tab on the screen [L2] in order to gain access to the following Menu:

8.3.1 Configuration

- Languages: EN: English; FR: French; DE: German; ES: Spanish; IT: Italian.
- Up/Down key: Programming the functions that can be activated from the head of the video camera.

The following can be selected: AWB/periph, AWB/Mirror, Brightness, Window, Details, Gamma, AWB/no action.

For the Brightness, Window, Details, and Gamma parameters, the up and down buttons make it possible to increase or decrease this parameter respectively.

- Blue shift: Shift in the white colorimetry towards the blue end of the spectrum.
- Red shift: Shift in the white colorimetry towards the red end of the spectrum.

8.3.2 Video

BRIGHTNESS	This parameter varies the general brightness of the image
WINDOW	Shutter window analysis: surface of the image on which the camera analyses the quantity of light received, in order to adjust the speed of the shutter to obtain a right exposure time. This function makes it possible to obtain a good exposure even when the image has zones with very different illumination (eg a dark orifice, and its very illuminated edges will require a small window).
DETAILS	Accentuating the level of detail of the image giving a sensation of relief. The ANTIMOIRE position settings have been optimized to apply an «anti-moire» filter when using a fibroscope
GAMMA	Modifies the dynamics of the dark areas: +1 reduces the difference between the light portions and the dark portions, therefore making these dark areas appear less so. -1 does the opposite: the dark portions are made even darker.

 The modified values are retained when the video camera is switched off.

8.3.3 Zero Reset Parameters

The video parameters for the specialities can be modified and are automatically saved.

In the event of a problem, this function initialises the video parameters using their default values.

9 Description of the back panel

9.1 Mains connectors/Fuses

The video camera is powered via the mains connector fitted to the back panel [C2], which must be connected to the mains using the cord provided with the camera. This connector features a fuses compartment with a cover [F1] as well as the master switch for powering up the device [S9].

When replacing a fuse, you must disconnect the video camera from the mains and use the same type of fuse. The "T" in "T800mA" stands for "Time-delay". Only use fuses with UL/CSA markings.

9.2 Video outputs

1. Standard:

This video camera features 3 video outputs: a composite output [C9], a Y/C output [C8], and an RGB output [C7].

A USB2 outlet [C5] which enables your video camera to be connected directly to a PC or to a video recorder. This will enable you to record video sequences of an operation using any video recording software.

 **Connect the video camera's USB2 output [C5] to a computer's USB2 input.**

 **Make sure that your PC or recorder complies with the standards in force (IEC 60601-1) and that it has been authorised for use in an operating unit.**

 **It is prohibited to carry out the operation with the USB2 output; a video output must be connected to the monitor.**

2. Option:

Two SDI outputs enable the video signal to be transmitted over long distances (to a remote screen, etc.).

9.3 Peripheral output

A peripheral output [C6] enables print-outs, screen shots, etc. via a 3.5 mm stereo jack cable, using the control wheel on the sensor. A «p» will be displayed on the screen when this function is activated from the head of the video camera.

 First of all, you have to configure the activation of this output via the Tool menu [cf. section 8.3].

 **Devices connected to the «VIDEO OUT» and «PERIPHERAL» sockets must comply with the IEC 60950 standard.**

9.4 Indications and symbols on the back panel

The indications and symbols on the back panel enable the video camera to be identified in accordance with the following international standards: IEC 60601-1, IEC 417 and ISO 980 [see section 15].

10 Procedure for disinfecting the sensor

Pre-disinfection cleaning	Disinfection	Serilisation
Immersion in an enzymatic cleaning solution (Aniosyme DD1, Hexanios G + R or surface cleaning ((Anios disposable cleaning cloth))).	Immersion in a glutaraldehyde solution (Anios Laboratory - Steranios 2%)	<ul style="list-style-type: none"> • Sterilisation using the Sterrad 100NX process standard cycle (47 minutes) • Sterilisation usign ethyle oxide

 Any other disinfection method is prohibited. Damage caused by these other methods cannot be covered by the manufacturer.

 Use of a sterile protective cover on the video camera is recommended regardless of the process adopted by the sterilisation service, with this applying throughout the surgical operation.

 The video camera cannot be placed in an autoclave.

 The video camera is not compatible with automatic washing and disinfecting machines.

 Existing alkaline solutions for the preliminary disinfection of certain medical devices are PROHIBITED for the purpose of preliminary disinfection of our video cameras.

 You absolutely must rinse copiously those parts that have come into contact with the disinfectant.

 In order to avoid scratching them, use non-woven compresses for drying the lenses.

 The procedures described in this section are provided as advice; under no circumstances may they take the place of official directives or recommendations.

 As decontamination depends on the products, methods and/or tools selected, this is entirely the responsibility of the staff in question.

11 After-sales service

No special maintenance on the video camera is required.

 **Reminder: Misuse is not covered by the guarantee.**

If a fault persists and it proves necessary to return the video camera to the after-sales service department, take care to send it in its original packaging.

Likewise, you should send back the whole video camera (control channel unit, sensor, objective, and cables). Please enclose a short explanatory note concerning the fault noted along with the dispatch note.

 **The equipment absolutely must be disinfected before it is returned for repairs.**

Upon return of the equipment, check its condition and state any reservations on the delivery slip if necessary, and then confirm them with the carrier via a registered letter within 48 hours. Should an equipment item sent by us become damaged during transport, the repair amount will be assigned to the carrier if reservations have been stated within this deadline, or to the recipient should this not be the case.

12 Troubleshooting

● 12.1 The stand-by indicator key [L1] does not light up when the device is powered up

- Check that the mains connector at the back of the video camera is properly connected to the mains and that the video camera's master switch is on («I»).
- Check that the fuses are in good condition (only use T800mA fuses - 250V - UL/CSA).

● 12.2 The stand-by indicator key [L1] light up and stay on, but no image appears on the screen

- Check that the sensor is properly connected to the control channel unit (if not, a bar test pattern will be displayed).
- Check that the video camera's control channel unit is properly connected to the monitor (cable in good condition, plugs properly inserted).
- If the video chain features several items, connect the video camera directly to the monitor in order to ascertain whether the issue is not something to do with one of the peripherals.
- Check that the monitor is powered up, that the right video input has been selected, and that the images settings are not in the minimum position (colour, light and contrast).
- Check for the presence of light by checking the light source, the light cable, and the optics (endoscope).

● 12.3 If the image is fuzzy

- Check that there is nothing fouling the lens of the video camera or endoscope.
- Check that the lens is properly adjusted.
- Check the "Flexible" or "Rigid" position in the Menu and select the one you want to use.

● 12.4 The image is extremely light or dark

- Check that the Brightness parameter is not set to the maximum or minimum setting.
- Check that the Window used is neither too large nor too small (in manual mode).
- Check that the light source is active.

13 Technical characteristics

- Class I device, BF type applied part
 - Sensor: 1 CCD 1/4"
 - Resolution: [752 x 585] PAL; [768 x 494] NTSC
 - Display: [625 lines/50 Hz] PAL; [525 lines/60 Hz] NTSC
 - Definition: >500 lines
 - Sensitivity: >1 lux
 - Signal/noise ratio: 52 dB (gain of 0 dB)
 - Electronic shutter: Automatic (1/50 at 1/100,000 sec.)
 - Browsing wheel
 - Programmable sensor keys functions
 - OLED screen
 - VDC compatible with PAL or NTSC sensor
 - Variable shutter (window) analysis surface
 - Contour accentuation
 - Anti-moire filter
 - Separate adjustment of the Blue and Red shift
 - Gamma adjustment
 - Menu language selection (French, Italian, Spanish, English, German)
 - Colour bar test pattern
 - White balance with automatic memory
 - Automatic saving of user settings parameters in the memory
 - 1 composite video output (1 VDC/75 Ω)
 - 1 Y/C video output ([1 VDC & 0.3 VDC/75 Ω]
 - 1 RGBS output (3 x 0.7 VDC & 0.3 VDC/75 Ω)
 - 1 USB2 output
 - 1 peripheral output
 - Power supply: (100-240 V~; 50 Hz/60 Hz)
 - Power consumption: 20 VA
 - 2 x T 800 mA fuses - 250 V UL/CSA
 - Equipotential socket
 - Dimensions of the control channel unit (w x d x h): 300 x 300 x 66 mm
 - Weight of the control channel unit: 1,500g
 - Operating temperature: +10°C to +40°C
 - Transport and storage temperature: -10°C to +45°C
 - Relative humidity: 20 to 85 %
 - Atmospheric pressure: 700 hPa to 1,060 hPa
 - Control channel unit not protected against splashing water (IPX0)
 - Not suitable for use in the vicinity of an inflammable anaesthetic mix with air, oxygen or nitrous oxide
 - Continuous service
 - Complies with European Directive 93/42/EEC
 - Complies with standards IEC60601-1, IEC 60601-1-2, IEC 60601-2-18, and EN 980
- Optional features:
- 2 SDI outputs

12, 17 AND 18MM LENSES

- C Mount
- Endoscope holder
- Weight: 74g
- Airtight

CHARACTERISTICS OF THE STRAIGHT SENSOR

- Sensor: 1 CCD 1/4"
- Front plug for the aluminium C Mount
- Airtight video camera head made from ABS
- Dimensions: L: 80mm; Ø: 30 mm
- Weight: 210g
- Cable length: 3 metres
- Airtight

14 Electromagnetic compatibility

Manufacturer's guide and declaration – electromagnetic emissions

The "CCD mono video camera" equipment is intended for use in the electromagnetic environment specified below.

The user will have to make sure that it is used in this environment.

Emission test	Compliance	Electromagnetic environment - Guide
RF emissions CISPR 11	Group 1	The "CCD mono video camera" reference equipment uses radio energy only for its internal functions. Consequently, its RF emissions are very low and are not likely to create interference with nearby electronic equipment. The "CCD mono video camera" reference equipment must be used in all locations other than household rooms and those directly connected to the low voltage power public distribution network used for supplying buildings for domestic purposes.
RF emissions CISPR 11	Class A	
Harmonic Emissions IEC 61000-3-2	Complies	
Voltage fluctuations / Flicker IEC 61000-3-3	Complies	

14 Electromagnetic compatibility

Manufacturer's guide and declaration – electromagnetic immunity

The "CCD mono video camera" equipment is intended for use in the electromagnetic environment specified below.

The user will have to make sure that it is used in this environment.

Immunity test	IEC 60601 Level of gravity	Compliance level	Electromagnetic environment - guide
Electrostatic discharges IEC 61000-4-2	± 6 kV in Contact ± 8 kV in air	± 6 kV ± 8 kV	The floor must be made from wood, concrete or tiling. If the floor is covered with a synthetic material, the relative humidity must be less than 30%.
Transient voltage in bursts IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input/output lines	± 2 kV ± 1 kV	The quality of the main power supply must be that of a typical hospital or commercial environment.
Voltage surges IEC 61000-4-5	Differential mode ± 1 kV Common mode ± 2 kV	± 1 kV ± 2 kV	The quality of the main power supply must be that of a typical hospital or commercial environment.
Dips, short interruptions and variations in the power voltage IEC 61000-4-11	<ul style="list-style-type: none"> • < 5 % UT – for 10 ms • 40 % UT – for 100 ms • 70 % UT – for 500 ms • < 5 % UT – for 5 s 	< 5 % UT 10 ms < 40 % UT 100 ms < 70 % UT 500 ms < 5 % UT 5 s	The quality of the main power supply must be that of a typical hospital or commercial environment. If the user of the "CCD mono video camera" reference equipment requires it to be able to continue to function during interruptions in the main power supply, it is recommended that the "CCD mono video camera" reference equipment is supplied by an uninterruptible power supply or a battery.
Magnetic field at the network's frequency [50/60 Hz]	3 A/m	± 6 kV via contact ± 8 kV via air	The magnetic field at the network's frequency must be at a level characteristic of a location in a typical hospital or commercial environment.

Note: UT is the rated value of the power supply voltage applied during the test.

14 Electromagnetic compatibility

Manufacturer's guide and declaration – electromagnetic immunity

The equipment «mono ccd camera» is designed to be used in the electromagnetic environment specified below. The user should make sure that it is used in this environment.

Immunity test	IEC 60601 Level of gravity	Compliance level	Electromagnetic environment - guide
RF conducted IEC 61000-4-6	3 Vrms 150 kHz at 80 MHz	3 V	<p>Mobile and portable RF communication equipment must not be used at a distance from the "CCD mono video camera" reference equipment (including the cables) that is less than the recommended separation distance, calculated utilising the formulas applicable depending on the transmitter's frequency.</p> <p>Recommended separation distance</p> $d = 1,16 \sqrt{P}$ $d = 1,16 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2,33 \sqrt{P} \text{ 800 MHz to 2.5 GHz}$ <p>where P is the maximum output power of the transmitter in Watts (W) assigned by the manufacturer of the transmitter and d is the recommended separation distance in metres (m). The field levels emitted by fixed RF transmitters, which are determined by an electromagnetic measurement of the site, must be lower than the compliance level in each frequency band. Disturbances may occur in the vicinity of devices featuring the following symbol: </p>
RF radiated IEC 61000-4-3	3 V/m 80 MHz at 2.5 GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency band applies.

Note 2: These recommendations may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection due to structures, objects and people.

a The field levels of fixed transmitters, such as the base stations of radio telephones (cell phones/cordless phones) and ground mobile radios, amateur radio, and TV, AM and FM radio communication, cannot accurately be assessed theoretically. To obtain the electromagnetic environment appropriate for fixed RF transmitters, site measurements must be carried out. If a field level measured in the operating environment for the "CCD mono video camera" reference equipment exceeds the applicable compliance levels above, then the proper functioning of the "CCD mono video camera" reference equipment must be checked. If abnormal functioning is noted, additional measures must be taken, such as reorienting or relocating the reference equipment.

b Above and beyond the 150 kHz to 80 MHz frequency band, the field level must be lower than 3 V/m.

14 Electromagnetic compatibility

Recommended separation distances between portable and mobile RF communication equipment and the «CCD mono video camera» reference equipment

The equipment «mono ccd camera» is designed to be used in an electromagnetic environment in which the radiated RF The “CCD mono video camera” reference equipment is intended for use in an electromagnetic environment where radiated RF disturbances are controlled. The user of the “CCD mono video camera” equipment can help prevent electromagnetic interference by maintaining a minimal distance between the mobile and portable RF communication equipment (transmitters) and the “CCD mono video camera” equipment, as recommended below, depending on the maximum output power of the communication equipment.

Maximum output power assigned of the transmitter W	Separation distance depending on the frequency of the transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1,16 \sqrt{P}$	$d = 1,16 \sqrt{P}$	$d = 2,33 \sqrt{P}$
0.01	0.116	0.116	0.233
0.1	0.366	0.366	0.736
1	1.16	1.16	2.33
10	3.66	3.66	7.36
100	11.6	11.6	23.3

For transmitters with a maximum output power rating that is not listed above, the recommended separation distance d in metres (m) can be determined by using the equation applicable to the transmitter’s frequency, where P is the maximum output power of the transmitter in Watts (W) assigned by the transmitter’s manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance provided in the higher frequency band applies.

Note 2: These recommendations may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection due to structures, objects and people.

15 Symbols

- S** Button
- L** Screen
- C** Socket
- I** Label (rear panel of the product)

SB Sensor button

Date of manufacture

Manufacturer

Protective grounding

Class 1 Class I product

Compliant with European directive 93/42/CEE

Warning

An equipotential grounding socket conductor other than protective grounding or a neutral conductor, for direct connection between the electrical equipment and the potential equalisation busbar. Kindly refer to standard IEC 60601-1, 3rd edition.

T UL/CSA time-delay fuses

Video output ports and peripheral control output port

Data input and peripheral output

BF-type device

Electronic and electrical equipment marketed after 13/08/2005.

This symbol indicates that this product is not to be treated with household waste.

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