

User manual  
Manuel utilisateur  
Manual del usuario  
Manuale utente  
Benutzerhandbuch  
Gebruikershandleiding  
Podręcznik użytkownika  
Instrucoes de funcionamento  
Bruksanvisning  
**Ръководство за употреба**  
**УПАТСТВО ЗА КОРИСТЕЊЕ**  
Uputstvo za upotrebu  
Příručka pro uživatele

**UBICAM**  
digital diagnosis set


INTUITIVE  
CONNECTION


**COMEg**  
medical technologies

Camera connector


● User manual	5
● Manuel utilisateur	25
● Manual del usuario	45
● Manuale utente	65
● Benutzerhandbuch	85
● Gebruikershandleiding	105
● Podręcznik użytkownika	125
● Instrucoes de funcionamento	145
● Bruksanvisning	165
● Ръководство за употреба	185
● УПАТСТВО ЗА КОРИСТЕЊЕ	205
● Uputstvo za upotrebu	225
● Příručka pro uživatele	244

# English

INTUITIVE  
CONNECTION

# Contents

1	INTRODUCTION	8
2	INTRODUCTION OF THE CAMERA	9
3	PRECAUTION FOR USE	10
4	REGULATORY RECOMMENDATION	11
5	INSTALLATION	12
6	USE ADJUSTEMENTS	14
7	RECOMMENDED PROCEDURE FOR DECONTAMINATION	16
8	AFTER SALE SERVICE	17
9	TECHNICAL SPECIFICATIONS	18
10	ELECTROMAGNETIC COMPATIBILITY	19

# 1 Introduction

Thank you for your confidence shown in us by your purchase of this device.

To take full advantage and comply with all necessary precautions, you should carefully read this manual and be well acquainted with its contents.

Sentences containing the symbol  are points requiring special attention.

Sentences containing the symbol  are for information.

To facilitate the installation and the use of the camera we wanted to make his handbook easier to use. Thus, references in the page of presentation in the back of the cover (in the form of [D1] for example) are proposed for you to visualize easily the parts of the product concerned.

This user manual is an integral part of the appliance. It must be made available to the user. Correct usage and proper hanling of the appliance is made possible by following this manual. The user bears the sole responsibility for damages that result from misuse.

# 2 Presentation of the camera

It is a waterproof, USB2 colour micro camera with detachable sensor, designed for consulting offices diagnosis directly from your PC computer. Its lightness, its auto-shutter, its performances in terms of sensitivity, resolution as well as its fidelity in the respect of colours make of it the ideal diagnosis tool.

It includes:

- A 1/4" CCD, C-Mount camera head, at the end of a 3m cable.
- A USB2 cable with its waterproofness cap.
- A CD including:
  - The patient management software
  - A User Manual.
  - A USB2 driver.

In option:

- A 12 or 17mm waterproof objective and its endoscope holder ring.
- A 5m USB2 extension cord.

**This equipment was delivered in a cardboard box. You must keep this element for further transport of the camera.**

 **For the United States, Federal law restricts the sale of this product to physicians or under their order.**

# 3 Precaution for use

## **Read the user manual.**

- Respect the use and storage conditions.
- Install the camera in a clean, dry and well ventilated place.
- The installation of the camera must be exclusively done by a competent technician trained by the manufacturer.
- To avoid risk of electrical shocks, fire, short-circuits or dangerous emanations, do not insert any metallic object in the apparatus.
- Do not expose the camera to water and do not store it in damp places without the waterproof cap on the USB connector [risk of electric shock].
- Use only the accessories supplied with the apparatus or proposed as option by the manufacturer.
- Do not use this camera near any flammable anaesthetics.
- This camera must not be used in a ionizing environment.
- This camera is not a sterile device.
- Stacking and location close to other equipment may negatively affect EMC performance.
- Never use corrosive or abrasive products to clean the device, but only the disinfectant agents as recommended in chapter 7.
- Do not expose the device to splash water or use in an environment with excessive humidity.
- Other cables and accessories may negatively affect EMC performance.
- Before every use, check the device for rough surfaces, sharp edges or protruding elements which could present a safety risk.
- The device is not explosion-proof. therefore, never use it together with flammable anaesthetics.
- Do not place heavy objects on the device.
- The devices that connect to the inputs/outputs must comply with IEC 60950-1.
- Prior to activation, ensure the apparatus is not damaged. Deteriorated cables and connectors, must be replaced immediately.
- It is advisable to have a second camera available in the operating room to use if an absence or degradation of performance is observed.
- Check the compatibility of your endoscope and your light source with the manufacturer before use.
- If other cables or accessories are used than the ones specified this may result in an increase emissions or a decrease of the noise immunity of the device.
- Never place the device near a heat source or in a location where it is exposed to vibration and/or shock.
- The surface temperature of the camera can reach 41°C (after a few minutes' use). Therefore, avoid any contact with this area.
- Not suitable for use in the presence of inflammable anaesthetic gas mixture with air, oxygen or nitrous oxide.
- In the area of installation, the presence of high interference fields (strong magnetic fields for example) should be avoided. These interference fields can result in malfunctions of the camera.
- This appliance must be used on individuals (patients) suitable for undergoing a procedure of endoscopy.
- The camera is intended for professional use in a consulting room / operating theatre.
- Do not drop the appliance. If the appliance is dropped, do not reconnect it but return it to your authorized distributor or directly to the after-sales department.
- Do not submit the device to excessive dust exposure.
- The device should only be opened by a qualified technician authorized by the manufacturer.
- This appliance may not be modified without authorization of the manufacturer. If the medical device is modified, an inspection and a test must be performed to ensure that the medical device complies with the safety instructions.

# 4 Regulatory recommendation

## 4.1 Conformity

This device has been designed and manufactured by a company with a certified quality system. It complies with the requirements of the European Directive 93/42/EEC related to medical devices. As a result, it complies in particular with the ad hoc electric safety (IEC) and Electromagnetic Compatibility (EMC) standards.

## 4.2 Electromagnetic interference and electrostatic discharges

Although this product complies with the EMC standards, it is possible under very special circumstances that it may interfere with other devices or be interfered by other devices or by an unfavourable electromagnetic environment.

This medical electrical equipment needs special precautions regarding EMC and needs to be installed according to EMC information.

Image transfers from the camera to a monitor may be susceptible to electromagnetic disturbances.

To prevent such situations, it is recommended **i**:

- To take care of the quality of the electric network (and particularly the grounding system of all the devices and carriages);
- To keep the device away from Electromagnetic sources (e.g., compressor, motor, transformer, HF generator, etc.).

## 4.3 Materio vigilance **i**

As for any medical system, this device is subject to the provisions pertaining to materio vigilance; any serious malfunction must be promptly communicated to the competent authorities and the manufacturer with the greatest possible precision.

For the address of the manufacturer please see the last page of manual.

## 4.4 End of life

This appliance is marked according to the European Directive 2002/96/EEC on Waste Electrical and Electronic Equipment (WEEE).

By ensuring this product is disposed of properly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for recycling of electrical and electronic equipment.

Disposal must be carried out in accordance with local environmental regulations for waste disposal.

For more detailed information about treatment, recovery and recycling of this product, please contact the local dealer to whom you purchased the product from and who will tell you how to proceed.

# 5 Installation

The purpose of this device is the image of the work area through an endoscope.

This medical device is intended to be used by a qualified surgeon during endoscopic explorations. Its installation requires no special training. Please refer to the instructions in this manual.

## 5.1 Installation of the camera on the computer

A minimum configuration of the computer is necessary.

Minimum System Requirements	Recommended System Requirements
Intel® Pentium iv – 1,3 ghz processor	Intel® core i5 processor
512 Mb of RAM	2 Gb of RAM
250 Gb of hard disk space	320 Gb of hard disk space
2 USB 2.0 Hi-Speed ports	4 USB 2.0 Hi-Speed ports
Operating system: Windows® xP Pro sP3	Operating system: Windows® 7 Pro sP1
Graphics Card: 32 MB of unshared video RAM with support for DirectX 9	Graphics Card with NVIDIA or ATI chipset / 512 MB of unshared video RAM with support for DirectX 9
USB chipset: Intel ou NEC® / Renesas®	USB chipset: Intel ou NEC® / Renesas®
Screen Resolution: 1024 x 768	Screen Resolution: 1280x 1024 or higher

A USB2 driver included on the CD-ROM is necessary to the functioning of the camera.

 **DO NOT PLUG THE CAMERA ON A USB PORT ON THE FRONT PANEL OF THE COMPUTER.**  
 **INSTALL THE SOFTWARE BEFORE PLUGGING THE CAMERA.**

## 5.2 Installation of the imaging software and of the drivers

Refer to the software installation manual located in the document directory of the Imaging software CD-ROM.

# 5 Installation

## 5.3 Optimizing the image

**i** Considering the disparity of video performances between computers, it is peremptory to adapt the adjustments of the system ( graphic board, screen,...) to obtain the best image quality.

This procedure must be performed once at the installation of the camera.

- > Plug the camera to the computers, connect the endoscope, the light cable and start the light source.
- > Launch the software with the shortcut  on the desktop
- > Press on the icon to start the video 
- > Launch an AWB: wrap the endoscope in a white gaze to avoid any dark area on the image
- > Shift the resolution of the video to 720 x 576 for a PAL camera [Go to Files, Configuration, camera, capture connection].
- > Adapt the parameters of : Brightness/Contrast/ Saturation/ Sharpness in the following manner:
  - Go to Files, Configuration, camera, filter of capture, Amp proc video
  - Adjust the contrast to have a bright image, without white saturation or noise.
  - Adjust the brightness to have a bright and noiseless image.
  - Adjust the saturation of the colors of the images to obtain realistic colors.
  - Adjust the sharpness to enhance details without increasing the noise.
- > If necessary, modify the adjustments of the monitor to adjust colors and brightness.
- > Save the parameters and leave.

**i** A mode «By default» is available to come back to factory settings.

These values are kept and reloaded at each starting of the system.

Before connecting your telescope or fibroscope on the camera, you must check its condition (clean distal lens, focus if necessary) and link it to your light source.

**Then check if the image produced with such a configuration is good. Power your light source and place the extremity of your telescope in the palm of your hand, thus simulating the basic conditions of endoscopy (beware of the working distance "filmed object-telescope" which must be similar to reality, i.e. some centimetres).**

1. Wait a few seconds in order to have a stabilised light;
2. Launch the Automatic White Balance (see following chapter);
3. Place the extremity of your telescope in the palm of your hand, thus simulating the basic conditions of endoscopy (beware of the working distance "filmed object-endoscope" which must be similar to reality, i.e. some centimetres). If you obtain a fuzzy image, modify the focus of the lens and if necessary of the endoscope. Check also if there is neither fog nor dust on your lens.

# 6 User adjustments

## 6.1 White balance

**i** Prior to performing an AWB, the camera must be plugged to the computer with the endoscope and the light source connected. If there is not enough light, the AWB will not be done.

- Film a white surface, a white gaze preferably and take care it is properly placed before the endoscope. (beware the produced image is clear), see Chap. 5.3.
- Press on the button of the white balance of the camera [white button] or of the software.
- A white square appears on the screen indicating a pending white balance, do not move until this square disappears, indicating the end of the white balance.

Each balance is memorised. Thus, the next time you turn on the camera [and if no parameters of your image chain were modified], you will find the same colours.

You need to do a new white balance if:

- You change the type of light source;
- You change your lamp (new one);
- You change the optics;
- You change the light cable;
- You are not satisfied with the reproduction of the colours.

**i** If you are not satisfied with the quality of the colours, think to check the adjustments of your monitor (Chroma, contrast,...)

## 6.2 The ICE function

With the new «ICE» (Intelligent Control Endoscope) device built in this camera, the white balance you just performed also allows to optimize all the settings whatever the endoscope used is.

For each white balance, this exclusive process analyses:

- the size of the endoscope to adapt the field of vision of the camera,
- the color temperature from the light source to determine its wavelength,
- the color temperature of the light going through the endoscope to restore realistic colors during the white balance.
- the minimum level of light to warn the user in case of failure.

# 6 User adjustments

## 6.3 Image capture

A short pressure on the black button will allow you to capture an image.

According to the parameter settings of your Imaging software, the image will remain frozen on the screen and you will have to press again on the same button to come back to live, or the image will automatically shift to the live mode.

In both cases, the captured image is automatically saved in the patient file.

## 6.4 Flexible/Rigid

**i** This function is used to guarantee an excellent image quality whatever the used endoscope.

**FLEXIBLE:** this function activates the ANTIMOIRE filter of the camera. It is used only for flexible endoscopes made with optic fibers.


**RIGID:** this function cancels the ANTIMOIRE filter of the camera and allows getting back the DETAILS of the image.

It is used for rigid and semi-rigid endoscopes.


A long pressure on the capture button (black) will activate the anti-moiré\*. A new long pressure will shift back to Rigid.

When the anti-moiré\* mode is activated, an F appears on the low right side of your screen.

\*Anti-moiré : correction of the natural default of aliasing due to the association of a CCD camera with an endoscope with optic fibres.

Moiré

Anti-moiré

# 7 Recommended procedure for decontamination

**⚠ CAUTION:** This procedure is only to be used if the cap of the USB connector is correctly closed.

**i** The camera sensor and the lens are designed for a cold decontamination by immersion (Steranios, Chlore Dioxyde, Hexanios, Tristel swipes....). All other methods of decontamination are prohibited. The damages due to these other methods could not be taken over by us.

**⚠ Only immerse the C-Mount sensor with its objective screwed on it.**

It is peremptory to respect the procedure consisting in cleaning the parts to immerse before decontamination and after decontamination, to rinse carefully with sterile water all the parts which have been in contact with the decontaminant.

**⚠ CAUTION !**

- It is peremptory to rinse abundantly the parts which have been in contact with the decontaminant.
- Use unwoven compresses for the drying of the distal lens so not to scratch it.
- The procedures described are advisory, they can, in no cases, be substituted to official recommendations or directives.
- The decontamination being linked to the products, methods and/or selected tools, remains entirely under the total responsibility of the people concerned.
- It is advisable to use a disposable, sterile, protective cover on the camera regardless of the decontamination procedure followed and throughout the surgical operation.
- The camera is not autoclavable.

## Disinfecting the lens surfaces

A dirty lens surface can interfere with observation. The lens surface must not have any marks or smudges.

To avoid scratching the lens surface, never use abrasive cloths or sponges to clean it.

To remove dust, dirt and other impurities not originating from the patient, wipe clean with a soft, lint-free cloth moistened with ethanol or 70% isopropyl alcohol. When cleaning, wear gloves for protection from chemical products.

Ensure that the equipment is completely dry before use.

**⚠ Use soft non-woven cloths to clean the lenses so as to avoid scratching them.**

# 8 After sale service

No special maintenance is necessary on the camera. Here is a list of common problems which are easy to troubleshoot.

For all other cases, contact our Maintenance Service that will help you in best delay. We remind you that the mishandlings are not guaranteed.

## 8.1 No image on screen

- Check if the camera is correctly connected and if the imaging software is opened and correctly installed (See Chap. 5).

## 8.2 Flickering image

- Check the minimal configuration of the computer in Chap.5.
- Keep only the imaging software active, close any other software.
- Do not plug the camera on the front panel of the computer.

## 8.3 The image is fuzzy uniformly white

- Check if there is no blur on the lens.
- Check if the sensor does not aim at a too lighted object.
- Check the focus of the lens.

## 8.4 The image is too clear or too dark

- Check if the « BRIGHTNESS» and «CONTRAST» parameters are well adjusted (See Chap.5).

**⚠**The product must be disinfected before returning it to repair.

If the problem remains and if you are led to send us back the camera, be careful to do it in its original packing. Likewise, you must send us the whole system (sensor, lenses, and cords). Be kind enough to join an explanatory note relative to the observed default to the expedition bill.

At the reception of your equipment, you will have to check its state and to enter a reservation on the delivery bill if necessary. You will then have 48 hours to confirm them by recommended letter addressed to the carrier. Beyond this delay, the carrier may deny these reservations. In the case where equipment shipped by us suffers damages during the travel, the amount of the repairing will be charged whether to the carrier if the reservation were entered in the delays, or to the consignee in the contrary.

Think to check rapidly the good functioning of the equipment which has just travelled.

Any other interventions being performed by our After Sales Service, no Service Manuals will be supplied.

# 9 Technical specifications

- BF Type camera
- Sensor: CCD 1/4" high sensitivity
- Resolution of the sensor: [752 x 582] PAL ; [768 x 494] NTSC
- Definition: 470 lines
- Sensitivity : 2 lux
- Signal/Noise Ratio : 52 dB
- Electronic shutter: Automatic [1/50 to 1/100 000]
- Variable surface of analyse of the shutter
- White Balance: Automatic
- Video enhancement «SHARPNESS»; electronic anti-moire filter for flexible fibrescope
- Modification of the video parameters
- Automatic detection of the colour temperature according to the used light source
- Automatic saving of the adjustments parameters
- 2 programmable buttons
- Cable length: 3 meters, Ø = 5 mm
- 1 USB 2.0 output
- Power supply: 5 V DC
- Power consumption : 2.5 VA
- Operating temperature range: +10°C / +35°C
- Operating relative humidity range: 30 to 75 %
- Transport and storage temperature range: -10°C / +45°C
- Transport and storage relative humidity range: 20% to 85%
- Operating, transport and storage atmospheric pressure range: 700hPa to 1060hPa
- Continuous Service
- Not adapted to a use close to a mixture of inflammable anaesthetic with air or oxygen or nitrous oxide.
- **CE** Comply with the European Directive 93/42/EEC, IEC 60601-2-18 and Amendment for CANADA and USA.

## FEATURES OF THE CAMERA HEAD:

- ABS waterproof camera head
- Front part of the C-Mount made of aluminium
- Size of the sensor: L = 87mm; Ø = 30 mm
- Weight of the sensor: 210 g

## ACCESSORIES:

- Objective F=12 or 17mm with endoscope holder
- Weight of the objective : 60g

## 1 CD-ROM including:

- The pilot of the camera
- The patient management imaging software
- Software tools.

# 10 Electromagnetic compatibility

## 10.1 Guide and declaration by manufacturer - Electromagnetic emissions

This camera is designed to be used within the electromagnetic environment specified below. The user should ensure that it is indeed used within this environment.

Emissions test	Compliance	Electromagnetic environment -Guide
RF emissions CISPR 11	Group 1	This camera only uses radioelectrical energy for its internal subsystems. Therefore, it emits very low RF energy and is not likely to interfere with nearby electronic devices.
RF emissions CISPR 11	Class B	
Harmonic emissions EN 61000-3-2	Not applicable	This camera can be used in all premises other than residential premises and premises connected directly to the public low voltage power distribution network used for powering residential buildings.
Voltage fluctuations / flicker EN 61000-3-3	Not applicable	

# 10 Electromagnetic compatibility

## 10.2 Guide and declaration by manufacturer - Electromagnetic immunity

This camera is designed to be used within the electromagnetic environment specified below. The user should ensure that it is indeed used within this environment.

Immunity test	CEI 60601 severity level	Compliance level	Electromagnetic environment - Guide
Electrostatic discharges EN 61000-4-2	± 6 kV via contact ± 8 kV via air	± 6 kV ± 8 kV	The floor should be in timber, concrete or tiling. If the floor is covered with a synthetic material, the relative humidity should be at least 30%.
Quick transients in bursts EN 61000-4-4	± 2 kV for power lines ± 1 kV for input/output lines	± 2 kV ± 1 kV	The quality of the main power supply should be that of a typical commercial or hospital environment.
Voltage shocks EN 61000-4-5	Differential mode ± 1 kV Common mode ± 2 kV	Not applicable	Not applicable
Brown-outs, brief power failures and voltage variations EN 61000-4-11	<5% UT – for 10 ms 40% UT – for 100 ms 70% UT – for 500 ms <5% UT – for 5 s	Not applicable	Not applicable
Magnetic field at network frequency (50/60 Hz)	3 A/m	3 A/m	The magnetic field at the network frequency should be at a level characteristic of a location within a typical commercial or hospital environment.
Conducted RF EN 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	Portable and mobile RF communication devices must not be used at a distance from this camera – including cables – that is less than the recommended separation distance calculated using applicable formulae according to the transmitter's frequency.  Recommended separation distance  $d = 1,16 \sqrt{P}$ $d = 1,16 \sqrt{P} \text{ 80 MHz à 800 MHz}$ $d = 2,33 \sqrt{P} \text{ 800 MHz à 2,5 GHz}$
Radiated RF EN 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	In which P is the maximum output power of the transmitter, in Watts (W), assigned by the manufacturer of the transmitter and the recommended separation distance (d) in meters (m). The field levels emitted by immobile RF transmitters – which are to be established by an electromagnetic measurement at the site – must be lower than the compliance level in each frequency band. Interference can occur nearby devices bearing the following symbol: 

NOTE : UT is the nominal value of the power voltage applied during the test.

NOTE 1 : At 80 MHz and 800 MHz, the higher frequency band applies.

NOTE 2 : These recommendations may not apply in all situations. The propagation of electromagnetic waves is modified by absorption and reflection due to structures, objects and people.

a The field level of immobile transmitters, such as radio telephone base stations (cellular and cordless) and mobile terrestrial radio systems, amateur radio systems, AM/FM radio communication systems, and TV systems cannot be evaluated theoretically with precision. To ascertain the electromagnetic environment due to immobile RF transmitters, a site measurement must be performed. If a field level measured within the environment in which this camera is used exceeds the above applicable compliance levels, check that this camera operates in a satisfactory manner. If abnormal operation is observed, additional measurements should be taken, such as reorientation or relocation of the reference system.

b Outside the frequency band of 150 kHz to 80 MHz, the field level should be less than 3 V/m

# 10 Electromagnetic compatibility

## 10.3 Separation distances recommended between portable and mobile RF communications systems and this camera

This camera is designed to be used within an electromagnetic environment in which radiated RF interference is controlled. The user of this camera can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication systems (transmitters) and this camera, such as recommended below, as a function of the maximum output power of the communications system.

Maximum assigned output power of the transmitter W	Separation distance as a function of the transmitter's frequency M		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1,16\sqrt{P}$	$d = 1,16\sqrt{P}$	$d = 2,33\sqrt{P}$
0.01	0.116	0.116	0.233
0.1	0.366	0.366	0.736
1	1.16	1.16	2.33
10	3.66	3.66	7.36
100	11.6	11.6	23.3

For transmitters of which the maximum output power is not listed above, the recommended separation distance d in meters [m] can be established by using the equation applicable for the transmitter frequency, in which P is the transmitter's maximum output power in Watts [W] assigned by the transmitter's manufacturer.

**NOTE 1:** AT 80 MHZ AND 800 MHZ, THE SEPARATION DISTANCE STATED FOR THE HIGHER FREQUENCY RANGE SHOULD BE APPLIED.

**NOTE 2:** THESE RECOMMENDATIONS MAY NOT APPLY IN ALL SITUATIONS. THE PROPAGATION OF ELECTROMAGNETIC WAVES IS MODIFIED BY ABSORPTION AND REFLECTION DUE TO STRUCTURES, OBJECTS AND PEOPLE.

## SYMBOLS


 Manufacturing date Manufacturer Comply with the European Directive 93/42/CEE Read the User Manual Dispositif de type BF Disposal of electric and electronic equipment marketed after 13/August/2005.  
This symbol indicates that the product cannot be treated with domestic waste.

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