


SOPRO 284

USER MANUAL	2
MANUEL D'UTILISATION	20
MANUAL DE USUARIO	38
MANUALE PER UTILIZZATORE	56
BETRIEBSANLEITUNG	74
GEBRUIKSHANDLEIDING	92

CONTENTS

1	FOREWORD	3
2	INTRODUCTION	4
3	SAFETY INSTRUCTIONS	5
4	STANDARDS AND REGULATIONS	7
5	INSTALLATION	8
6	OPERATION	9
7	DESCRIPTION OF BACK PANEL	10
8	RECOMMENDED DECONTAMINATION PROCEDURE	11
9	AFTER-SALES SERVICE & MAINTENANCE	12
10	TROUBLESHOOTING	13
11	TECHNICAL SPECIFICATIONS	14
12	ELECTROMAGNETIC COMPATIBILITY	15
13	KEY TO SYMBOLS	18

FOREWORD

Thank you for placing trust in us by purchasing this product.

To optimise the use of this device, it is essential to read this manual thoroughly and take all the necessary precautions.

Sentences bearing the  symbol refer to points which require special attention.

Sentences bearing the  symbol are for information purposes.

Our aim was to make the manual more practical to aid installation and use of the device. We have therefore provided references to the presentation page on the inside cover (for example, in the form of [D1]) to help you visualise the relevant parts of the product more easily.

PRESENTATION OF THE DEVICE

This is an LED light source intended for use by medical practitioners when performing surgical or diagnostic endoscopic procedures. It has an electronic iris with fibre optic sensing, an automatic thermal protection system a power level memory as well as an "iris" mode that allows the user to turn the light off.

Its ease of use, illumination power and precise manual adjustment of light intensity make this the perfect multidisciplinary medical instrument.

This device comprises:

- A power cord
- A user manual
- A quick start guide

Optional accessories:

- Universal light cable
- Adapters for universal light cable, source side (Storz and Olympus) and endoscope side (Storz, Olympus, Wolf, S&N, ACMI, Stryker)

We advise selecting the following diameter of light cable:

- Ø = 4.8 mm for endoscopes with diameter >4 mm
- Ø = 3.5 mm for endoscopes with diameter >=4 mm

⚠ For the USA and Canada, use the hospital-grade cord supplied, which must be plugged into a hospital-grade socket.

⚠ This equipment is delivered in cardboard packaging which should be kept for future possible transportation.

⚠ Under US Federal Law, sale of this product is restricted to medical practitioners or it may be sold on order of a medical practitioner.

SAFETY INSTRUCTIONS

Read the user manual.

- Comply with the conditions of use and storage conditions.
- The device may only be opened by a qualified technician authorised by the manufacturer.
- Do not insert any metal object into the device, in order to avoid any risk of electric shock, fire, short-circuit or hazardous discharge.
- Do not expose the device to an excessively humid environment or where it could be sprayed with water.
- Only use the accessories supplied with the device or optional accessories recommended by the manufacturer.
- Do not place any heavy objects on the device.
- This device is not sterile.
- If the power cord is damaged, immediately switch off the device. It is dangerous to operate this device with a damaged cord.
- To disconnect the cord, pull out the plug. Never pull on the cord itself.
- Disconnect the device from the mains if you do not expect to use it for a few days or longer.
- Ensure there is sufficient air circulating to prevent overheating inside the device: at least 15 cm all around the device. Do not cover the device and make sure the feet of the device are present.
- Do not use corrosive or abrasive products to clean the device, but only use the disinfectant liquids recommended in section 8.
- Before using the device, make sure it does not have any rough surface, sharp edges or protuberances that might cause safety problems.
- The surface temperature can reach 41°C (after a few minutes' use). Therefore avoid keeping this area in contact with the patient.

- To avoid any risk of electric shock, this device should only be connected to an earthed power supply.
- This device must not be modified without the manufacturer's authorisation. If the medical device is modified, an inspection and test must be performed to ensure the medical device complies with the safety instructions.
- Before use, check with the manufacturer that your endoscope is compatible with this device.
- This device must be used on individuals (patients) capable of undergoing an endoscopic procedure.
- Never look directly at the light aperture or the end of the light cable.
- Do not insert anything other than a light cable into the housing [C1] as this may damage the optical system.
- The use of cables or accessories other than those specified may result in increased emissions or decreased immunity of the device.
- After use of the source, when the fibre is being withdrawn from its light guide, the temperature on the metal coupling of the fibre is very high and can cause burns.
- The light power at the fibre outlet can cause eye damage. Take care when handling the fibre while the device is in use.
- Do not place the distal end of the light cable or the endoscope directly on the patient or any other inflammable material (surgical drapes, gauze, etc.) because it may be very hot and may cause burns.

 It is advisable to have a second light source available in the operating theatre, which can be used as a replacement if performance is seen to be reduced or absent.

 **This product is fitted with group 1 LEDs in accordance with standard IEC 62471. Do not look directly at the light to avoid any risk of eye damage.**

STANDARDS AND REGULATIONS

4.1 CONFORMITY

This product was designed and manufactured by a company which has a certified quality system.

It complies with the requirements of European Directive 93/42/EEC on medical devices.

It therefore conforms particularly to the standards of electrical safety (IEC) and electromagnetic compatibility (EMC) ad hoc.

4.2 ELECTROMAGNETIC INTERFERENCE AND ELECTROSTATIC DISCHARGE

Although this product conforms to EMC standards, under very particular circumstances it may interfere with other devices, be subject to interference from other devices or an adverse electromagnetic environment.

To avoid these situations it is advisable:

- to ensure the quality of the electrical network (especially the earthing of all appliances and trolleys),
- to distance the device from electromagnetic sources (e.g. compressor, motor, transformer, HF generator, etc.).

4.3 MEDICAL DEVICE VIGILANCE

Like any medical device, this product is subject to medical device vigilance provisions; any serious malfunction should therefore be the subject of a report to the competent authorities and to the manufacturer as soon as possible and as precisely as possible.

Manufacturer's contact details: please refer to the last page of the manual.

4.4 DEVICE END OF LIFE

This device bears the recycling symbol in compliance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

By correctly scrapping this device, you will help to prevent any harmful consequences for the environment and human health. The symbol  appearing on the device or on the accompanying documentation indicates that this product cannot under any circumstances be treated as household waste. Therefore, it should be given to a waste collection centre that handles recycling of electrical and electronic equipment. When scrapping this device, standards in force in the country of installation relating to waste disposal must be observed.

For more details about the treatment, recovery and recycling of this device, please contact your nearest dealer regarding the recommended procedure.

INSTALLATION

- This medical device is intended for use by a qualified surgeon during endoscopic exploratory procedures. Its installation does not require any special training. Please refer to the instructions contained in this manual.
- Place the device on a stable surface and pick up the various accessories that are essential to its operation.
- If you are using furniture, make sure this is adequately ventilated (at least 15 cm around the device).
- Connect the power cord to the mains socket on the device [C3].
- Set the switch [S3] to position 1, the device is then in standby mode (the green light [L1] should flash and a logo should appear in the display [L2]).
- Press the Standby button [S1] to switch on the device (the green light is lit and a logo indicating insertion of a light cable appears in the display [L2]). Pressing the Standby button again will reset the device to standby.
- Push the light cable into the aperture provided [C1]; you should hear a click as it engages. A logo indicating that the device is in "iris" mode appears.
- Connect the other end of the light cable to your endoscope.
Before setting up an endoscope (or fibroscope) in front of the camera lens, it is advisable to check the condition of the endoscope (distal lens or clean eyepiece).
- Press the adjustment knob [S2] or turn it in order to release the light.

 The light source is fitted as standard with an adapter for a Storz-type light cable. A version of this device that is compatible with Olympus-type light cables is also available.

 **If no light cable is connected to the device when the latter is switched on, a logo will display this information on the screen. Then you only need to connect up a light cable in order to see the "iris" logo and use the device.**

OPERATION

6.1 SETTING THE LIGHT INTENSITY

After insertion of the light cable, a logo indicating the "iris" mode will be displayed on the screen [L2].

Pressing or turning the adjustment knob [S2] enables the user to exit this mode and to release the light. A bar graph indicates the level of intensity will then appear [L2]

Turning the adjustment knob [S2] clockwise will increase the light intensity and turning it anti-clockwise will reduce the intensity.

An empty bar graph indicates that the source is completely dead while a full bar graph shows that it is delivering maximum light intensity.

 **The LED takes about 10 seconds to stabilise its colour temperature. It is advisable to wait for this length of time before setting the white balance on the camera.**

 **It is also advisable to do a white balance whenever there is a change of intensity setting.**

6.2 POWER LEVEL MEMORY

To put the LED source in "iris" mode, simply press the adjustment knob [S2]. Pressing the same button again will restore light power to the previous setting.

In the same way, if the LED source has been disconnected from its power supply it will start up again with the light power of the previous setting.

DESCRIPTION OF BACK PANEL

7.1 DESCRIPTION OF VENTS

The vents [A] located on the back panel of the device must not be obstructed as this can cause a risk of overheating.

The light source has an automatic safety device which stops the lighting function if the internal temperature is too high. A pictogram indicating whether this safety system is activated will be displayed on the screen [L2].

7.2 TRANSFORMER / MAINS SOCKET

The electrical supply to the device is via an IEC standard mains cord provided with the equipment.

7.3 MAINS SOCKET – FUSES

The electrical supply to the device is via the power socket located on the back panel [C3], which must be connected to the mains via the mains cord provided. This socket comprises a fuse flap [F] as well as the main 'on' switch [S3].

When replacing a fuse, it is essential to disconnect the product from the mains and use the same type of fuse. The T in 'T800mA' means 'time delay'. Only use UL/CSA labelled fuses.

7.4 INDICATORS AND SYMBOLS ON BACK PANEL

The indicators and symbols on the back panel of the product identify it as conforming to the international standards IEC 60601-1, IEC 60601-2-18, IEC 60417 and EN 980 (see section 13)..

RECOMMENDED DECONTAMINATION PROCEDURE

- ⚠ The LED source is a non-heat-resistant medical device which will not withstand immersion. It should therefore be disinfected using a non-woven cloth impregnated with a detergent-disinfectant product.
- ⚠ As decontamination is related to which products, methods and/or tools are chosen, it remains the sole responsibility of the personnel involved.
- ⚠ This device is not autoclavable.
- ⚠ The existing alkaline solutions for disinfection of certain medical devices are **PROHIBITED** for disinfecting this device.

AFTER-SALES SERVICE & MAINTENANCE

No special maintenance is required for this device.

⚠ Note: incorrect or faulty use is not covered by the guarantee.

If a problem persists and the product needs to be returned to our After-Sales Service, please be sure to send it to us in its original packaging.

Similarly, you should send us the whole device (control unit, cables). Please include a small note with your delivery note, giving details of the fault encountered.

⚠ It is essential to disinfect the equipment before return for repair.

When returning the equipment, check its condition and note any comments on the delivery note, if necessary. Please confirm these comments by registered letter to the delivery company within 48 hours. In a situation where equipment sent by us sustains damage in transit, the cost of repair will be invoiced either to the delivery company if the comments were sent within the deadline, or to the recipient if this has not been done.

In the event of any incident, please contact our Customer Service or our nearest sales representative.

TROUBLESHOOTING

10.1 THE LIGHT [L1] IS NOT LIT WHEN SWITCHED ON

- Check that the power cord is properly connected to the electricity supply and to the device and that the main switch [S3] on the back of the unit is in the 'on' position ("I").
- Check that the fuses [F] are in good condition (only use T800mA - 250V – UL fuses).

10.2 THE SOURCE IS LIT BUT THE FLOW OF LIGHT IS INADEQUATE

- Check the condition of your light cable and your optics.
- Check that the intensity setting is not too low.

10.3 THE FANS CONTINUE RUNNING AFTER THE LIGHT HAS GONE OUT

The light cable is disconnected.

The light source is fitted with a safety device which cuts the supply to the LED if the temperature inside the device is too high. A pictogram depicting the status of this safety system will be displayed on the screen. Once the temperature has fallen, the source returns to the Standby position and it can be re-used:

- Check that there is enough space around the source to ensure adequate cooling (15 cm all around the device).
- Check that nothing is obstructing the air vents on the back and underneath the source.

The light source is in "iris mode". Press the adjustment knob [S2] or turn it in order to release the light.

If the fault persists and the device needs to be returned to our After-Sales Service, please be sure to send it to us in its original packaging after the device has been disinfected.

10.4 LIGHT APPEARS / DISAPPEARS INTERMITTENTLY

Check that the light cable is fully engaged and the connector [C1] is properly screwed in place.

TECHNICAL SPECIFICATIONS

Main functions

- LED technology
- Colour temperature: 6500°K
- Typical LED life-time: 50,000 hours
- Compatible light cable: Storz-type or Olympus (optional)
- Electronic iris
- Light intensity adjustment: manual
- Anti-glare system with light cable sensor
- Standby mode
- Automatic thermal protection system
- Power level memory
- Supply voltage: 100-230VAC / 50 - 60 Hz
- Power consumption: 140 VA
- Dimensions: 300 x 66 x 300 mm
- Weight: 2400 g

Environment

- Operating temperature: +10°C / +40°C
- Operating humidity: 30% to 75%
- Storage and transport temperature: -10°C / +40°C
- Storage and transport humidity: 20% to 85%
- Operating, storage and transport atmospheric pressure: 800 hPa to 1060 hPa
- Not protected from vertically falling water (IPXO).
- Not suitable for use in the presence of an inflammable anaesthetic mixture with air, oxygen or nitrous oxide.

Environment

- Electrical safety class 1, type CF
- Complies with the European Directive 93/42/EEC - Class I
- Complies with the international standards IEC 60601-1; IEC 601-2-18; IEC 60417 and EN 980
- This product is fitted with group 1 LEDs in accordance with the standard IEC 62471. Do not look directly into the light to avoid any risk of eye damage.

ELECTROMAGNETIC COMPATIBILITY

12.1 GUIDE AND DECLARATION OF THE MANUFACTURER - ELECTROMAGNETIC EMISSIONS

The equipment known as "LED Source" is intended to be used in the electromagnetic environment specified below. Users should make sure it is used in this environment.

Emission trial	Compliance	Electromagnetic environment - Guide
RF emissions CISPR 11	Group 1	This product "LED Source" may be used in all domestic settings, including those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic current emissions EN 61000-3-2	Complies	
Voltage fluctuations / Flicker EN 61000-3-3	Complies	

12.2 GUIDE AND DECLARATION OF THE MANUFACTURER - ELECTROMAGNETIC IMMUNITY

The equipment known as "LED Source" was designed to be used in the electromagnetic environment specified below. Users must ensure that it is actually used in this environment.

Immunity trial	CEI 60601 severity level	Compliance level	Electromagnetic environment - Guide
Electrostatic discharges EN 61000-4-2	± 6 kV via contact ± 8 kV via the air	± 6 kV ± 8 kV	The floor should be wooden, concrete or tile. If the floor is covered with a synthetic material, the relative humidity should be at least 30%.
Fast transient bursts EN 61000-4-4	± 2 kV electric lines ± 1 kV linput/output lines	± 6 kV via contact ± 8 kV via air	The main power supply quality should be one of a traditional commercial or hospital environment.
Electrical shocks EN 61000-4-5	Differential mode ± 1 kV Common mode ± 2 kV	± 6 kV via contact ± 8 kV via air	The main power supply quality should be one of a traditional commercial or hospital environment.
Voltage dips, short interruptions and voltage variations EN 61000-4-11	<5% UT – for 10 ms 40% UT – for 100 ms 70% UT – for 500 ms <5% UT – for 5 s	± 6 kV via contact ± 8 kV via air	The main power supply quality should be one of a traditional commercial or hospital environment. If the user of this product requires it to continue operating during main power supply interruptions, it is recommended that this product be supplied by UPS or battery.
Magnetic field with the network frequency (50/60 Hz)	3 A/m	± 6 kV via contact ± 8 kV via air	The magnetic field with the network frequency should be at a characteristic level of a location (50 / 60 Hz) in a traditional commercial or hospital environment.

NOTE: *Ur* is the power voltage nominal value applied during the trial.

GUIDE AND DECLARATION OF THE MANUFACTURER - ELECTROMAGNETIC EMISSIONS

The equipment known as "LED Source" is intended to be used in the electromagnetic environment specified below. Users should make sure it is used in this environment.

Safety trial	CEI 60601 severity level	Compliance level	Electromagnetic environment - Guide
Conducted RF EN 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	<p>Portable and mobile RF communication equipment should not be used at a distance from this product, including cables, that is less than the recommended separation distance, calculated by applying the formula that corresponds to the transmitter frequency.</p> $d = 1,16 \sqrt{P}$ $d = 1,16 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,33 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$
Radiated RF EN 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	<p>where P is the maximum output of the transmitter in watts (W) rated by the transmitter manufacturer and (d) is the recommended separation distance in metres (m). The field levels emitted by the fixed RF transmitters – which must be determined by an electromagnetic measurement in situ – should be lower than the compliance level in each frequency band. Interference may occur in the vicinity of devices bearing the following symbol: </p>

NOTE 1 : At 80 MHz and 800 MHz, the higher frequency band should be used.

NOTE 2 : These recommendations may not apply in every situation. Electromagnetic wave propagation is modified by the absorption and reflection due to structures, objects and persons.

For transmitters whose maximum power output is not listed above, the recommended separation distance d in metres (m) can be determined by using the equation applicable to the transmitter frequency, where P is the maximum power output of the transmitter in watts (W) as rated by the transmitter manufacturer.

12.3 RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATION SYSTEMS AND THIS PRODUCT

The equipment referred to as “LED Source” was designed to be used in an electromagnetic environment in which the RF interference emitted is controlled. Users of this device can help to avoid electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and this device, as recommended below.

Rated maximum power output of transmitter W	Maximum separation distance depending on transmitter frequency M		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1,16 \sqrt{P}$	$d = 1,16 \sqrt{P}$	$d = 2,33 \sqrt{P}$
0.01	0.116	0.116	0.233
0.1	0.366	0.366	0.736
1	1.16	1.16	2.33
10	3.66	3.66	7.36
100	11.6	11.6	23.3

NOTE 1 : At 80 MHz and 800 MHz, the higher frequency band should be used.

NOTE 2 : These recommendations may not apply in every situation. Electromagnetic wave propagation is modified by the absorption and reflection due to structures, objects and persons.

For transmitters whose maximum power output is not listed above, the recommended separation distance d in metres (m) can be determined by using the equation applicable to the transmitter frequency, where P is the maximum power output of the transmitter in watts (W) as rated by the transmitter manufacturer.

KEY TO SYMBOLS

S Button

L Light

C Connection / socket

I Label (on the product)

 Date of manufacture

 Manufacturer

Class 1 Class I product

 Complies with European Directive 93/42/EEC

 Warning

 Symbol for potential equalisation

 **T** Time delay fuses UL/CSA

 CF type device

 Electronic and electrical equipment marketed after 13/08/2005.
This symbol indicates that this product must not be treated with household waste.

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