SOPRO
640



USER MANUAL
MANUEL D'UTILISATION
MANUAL DE USUARIO
BETRIEBSANLEITUNG
MANUALE PER UTILIZZATORE

SOPRO
COMEG

SOPRO 640-45L

SOPRO 640-30L

SOPRO 640-20L

External filter

USER MANUAL	P.2
MANUEL D'UTILISATION	P.20
MANUAL DE USUARIO	P.38
BENUTZERHANDBUCH	P.56
MANUALE DELL'UTENTE	P.74

ENGLISH

CONTENTS

1	INTRODUCTION3
2	SAFETY INSTRUCTIONS3
A.	USE	
B.	CONTRAINDICATION	
C.	WARNING	
3	STATUTORY RECOMMENDATIONS4
A.	COMPLIANCE	
B.	ELECTROMAGNETIC COMPATIBILITY	
C.	MEDICAL DEVICE VIGILANCE	
D.	END OF LIFE	
4	DESCRIPTION OF THE FRONT PANEL AND BACK PANELS5
5	SPECIAL FEATURES6
A.	CONTROLS FROM THE VIDEO CAMERA HEAD (S397 ONLY)	
B.	AUTONOMY FUNCTION WITH GAS CYLINDER	
C.	GREATER SAFETY OF USE	
D.	LOW FLOW	
E.	HIGH FLOW	
F.	AUTOMATIC LOW/HIGH FLOW SWITCHING	
G.	EXTERNAL PRESSURE RELIEF VALVE	
H.	MONITORING OF THE CYLINDER PRESSURE	
I.	MONITORING OF THE CENTRAL GAS WALL CO ₂ CONNECTION	
J.	LCD DISPLAY	
K.	CENTRAL CO ₂ GAS WALL SUPPLY AND CO ₂ BOTTLE CAPABILITY ON THE SAME DEVICE	
L.	PURGE FUNCTION	
M.	OBSTRUCTED TUBING MESSAGES	
N.	PREVENTIVE MAINTENANCE MESSAGES	
6	INSPECTION AND INSTALLATION IN THE OPERATING THEATRE8
A.	CHECKING THE DEVICE	
B.	ANCILLARY EQUIPMENT	
C.	OPTIONAL	
D.	ELECTRICAL CONNECTION	
E.	EQUIPOTENTIALITY	
F.	CONNECTING THE EXTERNAL FILTER	
G.	CONNECTING THE GAS CYLINDER	
H.	CONNECTING TO CENTRAL GAS WALL CO ₂ CONNECTION	
I.	DISCONNECTING THE GAS CYLINDER	
J.	DISCONNECTING FROM CENTRAL GAS WALL CO ₂ CONNECTION	
K.	CONNECTING THE COMMUNICATION CABLE	
7	CREATING AND MAINTAINING PNEUMOPERITONEUM9
A.	STARTING THE INSUFFLATOR	
B.	CONNECTING - DISCONNECTING THE PATIENT TUBING	
C.	MODIFYING PARAMETERS VIA THE MENU	
D.	PRESSURE INITIAL SETTING	
E.	STARTING INSUFFLATION	
F.	CONTROL	
G.	OVER-PRESSURIZATION	
H.	MONITORING THE FLOW	
I.	MONITORING THE VOLUME OF GAS USED	
J.	MONITORING THE VOLUME OF GAS AVAILABLE	
K.	END OF OPERATION	
8	FAULTS AND WARNING SIGNALS14
9	CLEANING / SERVICING / MAINTENANCE15
A.	CLEANING / DECONTAMINATION	
B.	SERVICING / MAINTENANCE	
10	SYMBOLS USED15
11	TECHNICAL SPECIFICATIONS16
12	ELECTROMAGNETIC COMPATIBILITY17

INTRODUCTION

Thank you for the confidence you have shown in buying this device.

Sections marked with the symbol  deal with points that require particular attention.

Sections marked with the symbol  are just for information.

SAFETY INSTRUCTIONS

 Please read this manual carefully before using the device, pay particular attention to the sections dealing with safety conditions and environmental conditions related to the device.

A. USE

This device is designed to create pneumoperitoneum (distension of the abdomen with CO₂ gas) as part of a diagnostic and/or operative laparoscopy procedure. Please refer to your laparoscope user manual for information on specific use.

The aim of this user manual is not to tell the user how to perform a laparoscopy but simply to explain how to use this insufflator.

FEDERAL LAW:

Federal Law in the United States restricts the use of this device solely to, or under the supervision of, a surgeon.

B. CONTRAINDICATION

 Use of this device in cases of intra-abdominal distension is contraindicated at all times when laparoscopy is contraindicated. Please refer to your laparoscope user manual for absolute and relative contraindications.

 This device is contraindicated for hysteroscopic insufflation; on no account it must be used in cases of intrauterine distension.

Note: the distension pressure of an insufflator for laparoscopy must not exceed 25 mmHg

C. WARNING**• Metabolic acidosis and resulting cardiac irregularity**

Avoid prolonged intra-abdominal pressures above 20 mm of mercury. This would cause one of the following risks:

- Reduction of respiration with compromised diaphragmatic movement
- Reduction of venous return
- Reduction of cardiac output
- Acidosis

Excessive absorption of CO₂ is due to either a too high flow rate or excessive pressure, or both. The abdomen can be sufficiently distended by a pressure between 10 and 15 mm of mercury. It is rarely necessary to select an abdominal pressure higher than 15 mm of mercury. At these levels, the extent of intravascular penetration should be low. Pressures higher than 20 mm of mercury are hardly ever necessary; they will increase the quantity and the speed of intravascular penetration. Adequate respiration helps to prevent problems associated with CO₂.

• Idiosyncratic reactions

For patients suffering from microdrepanocytic diseases or pulmonary insufficiency, the use of these devices can present an increased risk of metabolic imbalance related to the excessive absorption of CO₂.

• Hypothermia

High-flow insufflators involve a potential risk of hypothermia; therefore we recommend using a heating system to keep the patient's temperature stable.

• Gas flow

Surgical procedures must be carried out with insufflators able to reach flows between 4 to 10 l/min. Insufflators providing lower flow must only be used for diagnostics procedures.

• Gas coagulation system

Some medical devices delivering gas to the interior of the peritoneal cavity during surgical laparoscopies (gas lasers, for example) may lead to abdominal over-pressurization. If venous sinuses are formed, triggering an embolism, the use of a secondary gas source that is far less soluble in the blood than CO₂ can cause a rapid rise in intra-abdominal pressure. We recommend avoiding the use of these coagulation systems during laparoscopic procedures, as the manufacturer could not be held responsible for any incident associated with their use.

- **Bacterial filter**

The use of a hydrophobic bacterial filter is essential for preventing cross-contamination with the patient

- **Please follow the conditions of use and storage.**

- **Do not expose the device to too much dust.**

- **This device must only be opened by a technician approved by the manufacturer.**

- In order to prevent any danger of electrocution, fire, short circuits or dangerous emissions, do not insert metal objects into the device.

- Do not expose the device to water spray or splashes, and do not install it in a very humid area.

- Use only the ancillary equipment delivered with the device or suggested by the manufacturer as an option.

- Do not connect this device to any other device.

- This insufflator must not be used in the presence of flammable anesthetics.

- This device has not been designed for use in an ionizing environment.

STATUTORY RECOMMENDATIONS

A. COMPLIANCE

This device was designed and manufactured by a company with a certified quality system.

It complies with the requirements of European directive 93/42/EEC on medical devices.

Consequently it complies in particular with the specific standards relating to electrical safety (IEC) and electromagnetic compatibility (EMC).

B. ELECTROMAGNETIC COMPATIBILITY

Although this product complies with standards on electromagnetic compatibility (Appendix XII) it is possible that, in very particular circumstances, it could interfere with other devices, or even that other devices or an adverse electromagnetic environment could interfere with it.

To prevent similar situations, we recommend to:

- check the quality of the electrical network (and in particular the earthing of all the devices and trolleys)
- keep the device away from electromagnetic sources (for example mobile phones, GSM phones, compressors, motors, transformers, HF generators in general etc).

C. MEDICAL DEVICE VIGILANCE

Like any medical device, this device is subject to the provisions related to medical device vigilance; therefore any serious malfunction must without fail be promptly reported to the proper authority and to the manufacturer, giving as much detailed information as possible.

Please refer to the last page for the manufacturer's contact details.

D. END OF LIFE

This device carries the recycling symbol in accordance with European directive 2002/96/CE on Waste Electrical and Electronic Equipment (DEEE or WEEE).

By following the correct procedures for disposing of this device you will be contributing to the prevention of any consequences that could be harmful to the environment or human health.

The symbol  to be found on the device and/or in the accompanying documentation indicates that this product may on no account be treated as household waste. It must therefore be delivered to a waste collection point for the recycling of electrical and electronic equipment.

For rejecting the device, please observe the standards applied in the country of installation.

To obtain complete details about the treatment, recovery and recycling of this appliance, please contact your closest dealer to know how to proceed.

DESCRIPTION OF THE FRONT PANEL AND BACK PANELS

Reference on front panel	Description
A	5.7-inch LCD screen
B	Warning light
C	Tubing connector for patient insufflation
D	Tubing connector for external desufflation
E	Low/high flow selector button
F	Start/stop insufflation
G	Selection and validation thumbwheel
H	Start/standby button
I	Start/standby LED

Reference on back panel	Description
1	CO ₂ Input
2	On/Off switch
3	Fuse housing
4	Mains socket
5	Equipotential socket
6	Labeling
7	Connector for remote control and services

SPECIAL FEATURES

A. CONTROLS FROM THE VIDEO CAMERA HEAD (S397 ONLY)

From one of the three programmable buttons located on the camera head, it is possible to control the following functions of the S640 insufflator:

START (in low flow mode), HGH FLOW, STOP

A special communication cable must be connected between the S397 and the S640 (See S397 user manual for more information).

B. AUTONOMY FUNCTION WITH GAS CYLINDER (NOT AVAILABLE ON ALL MODELS)

This patented function gives the Remaining Time Autonomy (RTA) of the CO₂ cylinder. It helps the surgery team to better manage:

- the surgery timing,
- the bottle exchange.

The autonomy appears when the CO₂ cylinder pressure will be below 33 bar and when the autonomy will be equal or below 99 minutes.

When autonomy will be equal to 20 minutes 4 sound alarms will warn the surgery team and thus every 5 minutes until the end.

Autonomy can remain at a fixed value during a short or a long period, it will depend on many parameters.

C. GREATER SAFETY OF USE

• Pre-heating system

Heats the gas inside the device.

• Autotest

Automatic calibration of the device and control of its basic components in less than one second.

• Detection of tubing

The device will only start if a tube is both connected to outlet and external pressure relief valve (for safety reasons it is impossible to connect a single tube)

• Automatic regulation of the flow rate

This insufflator automatically regulates the flow rate to maintain an abdominal pressure equal to the pressure setting.

• An additional pressure sensor

If the measurement circuit fails, consistency of the measurements is continuously monitored; this means that the insufflation cycles can be interrupted if there is the slightest cause of doubt.

• A high-pressure relief valve

If over-pressurization occurs at the level of the high pressure regulator, a safety valve is available to limit any risk.

• Bacterial filter

Delivered inside the disposable tubing, this filter prevents risks from the insufflator to the patient

D. LOW FLOW

In this mode, the insufflation flow is limited to 2 l/min so as to create the pneumoperitoneum. This flow is not sufficient to regulate intracavitory pressure in the event of major leakages, switch in this case to high flow mode.

By default, the insufflator starts in low-flow manual mode. The surgeon can choose the automatic or the manual mode to activate the high flow.

E. HIGH FLOW

When the pneumoperitoneum is created, activating this mode permits a maximum flow between 20 and 45l, depending on the model. This capacity then makes it possible to compensate for any type of leakages.

F. AUTOMATIC LOW/HIGH FLOW SWITCHING

If automatic switching is selected from the menu, the device starts insufflating in low-flow mode, then switches to the high-flow position once blow-off pressure is reached.

G. EXTERNAL PRESSURE RELIEF VALVE

Our insufflators have a valve that allows CO₂ overpressure outside the device in order to avoid:

- Over- pressurization in the pneumoperitoneum if the pressure is more than 2 mmHg greater than the pressure set for more than 1 second and/ or if the pressure is more than 5 mmHg greater than the pressure set for more than 3 seconds.
- any return of fluid to the device's interior, thus preventing any risk of deterioration or contamination.

H. MONITORING OF THE CYLINDER PRESSURE

This last generation of insufflators is equipped with a system for measuring the cylinder pressure displayed in continuous.

Moreover, when the cylinder pressure reaches 40 bars a updated cylinder icon is displayed,  is flickering.

As soon as the cylinder pressure reaches 20 bars, the insufflator informs the user via:

- A message on the display in the Alarm area "End of bottle" and a sound alarm.
- A yellow light displayed

 From, 10 bars insufflator is not allowed to re-start

 When AUTONOMY function is ordered, cylinder pressure messages in the Alarm area are replaced by AUTONOMY messages.

 Use only with medical grade CO₂ for coelioscopy complying with European standards.

I. MONITORING OF THE CENTRAL CO₂ GAS WALL PRESSURE (NOT AVAILABLE ON ALL MODELS)

This last generation of insufflators is equipped with a system for measuring the central gas wall CO₂ pressure displayed in continuous.

As soon as the central gas wall CO₂ pressure reaches 2 bars, the insufflators inform the user via :

- A message on the display in the warning area "input pressure".
- A yellow light displayed.

 From, 2 bars insufflator is not allowed to re-start.

 From, 10 bars insufflators will automatically switch in "bottle" mode for safety reasons.

 Use only with central GAS WALL CO₂ supply complying with European standards.

J. LCD DISPLAY

In order to communicate better with the user and to improve the ease of use, we have provided our latest generation of insufflators with a 5.7-inch LCD display. This allows the following to be displayed in 5 different languages:

- warnings and error messages
- interactive menu allowing the device to be customized
- indication of empty cylinder
- indication of pressure set in mm Hg
- indication of intracavitory pressure in mm Hg
- indication of flow in l/min
- indication of volume of CO₂ used in liters
- indication of status: in start/stop and "high flow"/"low flow"
- CO₂ bottle on central wall supply icon

K. CENTRAL CO₂ GAS WALL SUPPLY AND CO₂ BOTTLE CAPABILITY ON THE SAME DEVICE (NOT AVAILABLE ON ALL MODELS)

This device is able to be connected to CO₂ network from 3.5 to 5 bars as well as CO₂ bottles. Inside the MENU simply select "WALL" or "BOTTLE" inside "CO₂ SOURCE" and it will be memorised

If a CO₂ bottle is connected when "WALL CO₂ source" is activated the device will automatically switch in "BOTTLE" mode for safety reasons. (Bottle Pressure > 10 bars)

L. PURGE FUNCTION

When a CO₂ bottle is connected there is the possibility across the MENU to purge the high pressure hose before disconnecting it from the bottle.

M. OBSTRUCTED TUBING MESSAGES

A message ears in the Alarm area when the S640 detects an obstructed tubing.

N. PREVENTIVE MAINTENANCE MESSAGES

When the device is powered on and during few seconds a message: "2000 surgeries before maintenance" is displayed, this number will decrease time after time and finally the message "Maintenance required" will be displayed.

 The maintenance will be made before 2000 surgeries or before one year of use it depends of the first event.

INSPECTION AND INSTALLATION IN THE OPERATING THEATRE

A. CHECKING THE DEVICE

Any deterioration, immediate malfunction or missing ancillary equipment must be immediately reported or confirmed by the carrier and/or dealer by registered letter so that it can be covered by insurance if necessary.

Please keep the original device packaging for its yearly maintenance return. In this event, please enclose a document to include your name, address and the reason for returning the device (problem encountered)

⚠ Place the insufflator on a flat and dry surface. This insufflator must be horizontal when used and must be placed higher than the patient to prevent any liquid from flowing back through the tubes.

B. ANCILLARY EQUIPMENT

Your insufflator is delivered with the following ancillary equipment:

1	user manual
1	main power cable 2.50 m long
1	open-ended spanner for the high-pressure hose connector and external filter connection
1	external 15µm CO ₂ filtre

C. OPTIONAL

- Pin connector, type DIN 477-1 no. 6 or US 7/16"
- High-pressure flexible CO₂ feed hose (Pmax: 200 bar), 1m long.
- Central gas wall supply low pressure CO₂ hose, 7m.
- Communication cable (to connect between S397 and S640)

D. ELECTRICAL CONNECTION

⚠ Check that the features of the main power supply include 100 to 230VAC and 50 Hz - 60 Hz with a required power of 75VA.

Connect the power cable on the rear panel in [4].

Connect the power cable to the power supply.

⚠ In the United States and Canada the "Hospital Grade" power cable supplied should be used, and it must, without fail, be connected to the "Hospital Grade" mains connection.

⚠ This insufflator is protected by two 2.5 amp time-lag supply fuses (UR) on the rear panel in [3]; do not use UR fuses of different sizes that are not certified.

⚠ To avoid any electric shock risk, the device must be connected only to the main network equipped with earth protection.

E. EQUIPOTENTIALITY

The safeguards against discharge of static electricity (ESD) are as follows:

- application of equipotential bonding on all devices is to be connected; connect the voltage compensator cable to the connector [5] on the rear panel.
- exclusive use of the ancillary equipment mentioned

Applications of equipotential bonding on all devices are also to prevent contact with a dangerous voltage if electrical insulation fails

F CONNECTING THE EXTERNAL FILTER

- Screw manually the external filter on the connector [1] on the rear panel.
- Tights it with the delivered wrench.

G. CONNECTING THE GAS CYLINDER

- ⚠ Maximum permitted operating pressure (60 bar)**
- ⚠ Never begin a surgical operation without being sure that there is a spare CO₂ cylinder available. Use only medical grade CO₂ for celioscopy complying with European standards.**
- Fit the CO₂ gas cylinder in an upright position so that it is firmly in place with its top pointing upwards, if possible on the same trolley as the insufflator.
- ⚠ Never use the insufflator if it is connected to a cylinder placed horizontally or one which has its top pointing downwards.**
- Ensure that a seal on the high pressure cable on the external filtre of the cylinder side is present (DIN only)
- Use the open-ended spanner provided to tighten the connector linked to the cylinder outlet, then the one linked to the external filter of the insufflator inlet port [1].

H. CONNECTING TO THE CENTRAL GAS WALL CO₂

Connect manually the low pressure hose to the wall connector.

- ⚠ Minimum permitted operating pressure 2 bars.**
- ⚠ Use only with central GAS WALL CO₂ supply complying with European standards.**

I. DISCONNECTING THE GAS CYLINDER

- ⚠ Before unscrewing the high-pressure hose:**
 - Check that the cylinder valve is closed, and then lower the pressure inside the system by gently loosening it.
 - Lower the pressure inside the system by gently loosening it.
- OR
- Use the Purge function (inside MENU, select PURGE and follow the screen messages) and unlock the hose.

J. DISCONNECTING FROM CENTRAL GAS WALL CO₂ CONNECTION

- Disconnect manually the low pressure hose from the wall connector

K. CONNECTING THE COMMUNICATION CABLE

Turn off the devices.

Screw the 3 pins connector of the communication cable on the rear side of the S640 [7]

Screw the 7 pins connector of the communication cable on the rear side of the S397

Turn on the devices

CREATING AND MAINTAINING PNEUMOPERITONEUM

⚠ This insufflator is intended only for diagnostic or operative laparoscopy. Any use outside this area constitutes misuse of the product for which the user will be responsible and for which the manufacturer declines all responsibility.

This device must be used only by qualified staff. The surgeon and the anesthetist will be permanently responsible for the device, with the anesthetist guaranteeing particular supervision of blood gases. The safety features with which this device is fit will not in any way relieve the medical staff of their responsibility for the permanent monitoring and supervision of the patient.

A. STARTING THE INSUFFLATOR

User action	Device response
Operate the switch [I] on the rear panel [2]	The device starts to run, the indicator light [I] on the front panel flashes
Open the gas cylinder	
Press the button [H] on the front panel.	The device runs its AUTO TEST The indicator light [I] on the front lights up. The welcome page appears indicating the type of insufflator Followed by the menu

[i] Remaining Interventions message appears for preventive maintenance.

[i] The AUTOMATIC TEST should run if tubing is not connected. It only takes a minute.

You may now connect the tubing as described in the following chapter

B. CONNECTING - DISCONNECTING THE PATIENT TUBING

⚠ The manufacturer can not be held responsible for any reaction or dysfunction associated with the use of damaged or unsuitable tubing.

⚠ Do not use tubing if the packaging is damaged.

⚠ Sterile tubes are disposable; do not re-sterilize.

⚠ Use of a hydrophobic bacterial filter is essential for the prevention of cross-contamination with the patient. Replace it for each patient.

Open the first set of packaging film on the tubing. Hand the second set to sterile personnel.

⚠ Keep the patient's sterile tube in the operating theatre.

Open the second sterile pouch and connect the tubing to the insufflator's gas outlet.

Fit the other piece of tubing into the external desufflation actuator

⚠ The tube needs to be freely inserted without kinks, and must not be obstructed.

- Connect the Luer lock to the trocar port on the patient side.

C. MODIFYING PARAMETERS VIA THE MENU

User action	Device response
<p>At the start of the procedure use the thumbwheel to select the MENU and press it to enter in the MENU.</p> <p>To modify a parameter, use the thumbwheel to select it [G].</p> <p>Press the thumbwheel to select:</p> <ul style="list-style-type: none"> • a language from the following: English, French, German, Italian and Spanish. The menu's default language is English. • Automatic or manual low/high flow switching • reset the volume of CO₂ used to zero. • reset the chronometer to zero • contrast control • purge function • CO₂ supply (bottle or wall) <p>Turn the thumbwheel to change the values and press the thumbwheel to validate.</p> <p>To exit, select "return" and validate.</p>	<p>The parameter is highlighted in the display</p>

⚠ Automatic switching must be selected by order of the surgeon.

D. PRESSURE SETTING

⚠ This operation must be carried out by or under the control of a surgeon.

Use the thumbwheel to select the pressure intra abdominal press it to begin your choice and turn it to select the required intra-abdominal pressure between 0 and 25 mmHg.

[i] The pressure commonly used is 12 mmHg.

[i] You will be asked as a safety precaution to confirm pressure set above 15 mmHg.

E. STARTING INSUFFLATION

Insufflation will start once the button [F] on the front panel is pressed. RUN appears on the screen when the insufflation is running.

[i] The MENU can not be accessed while the device is operating; to access it, insufflation has to be stopped by using the button [F].

Insufflation starts in low mode, 2 l/min, so that pneumoperitoneum is created under the safest conditions.

When the low-flow is activated : LOW appears in larger on the screen.

High-flow mode is started manually by pressing the button [E] on the front panel, or automatically if the mode HIGH is activated (in the MENU).

When the high-flow is activated : HIGH appears in larger on the screen.

F. CONTROL

Once the cavity is created, the insufflator keeps pneumoperitoneum at the pressure selected and corrects CO₂ leakage immediately.

⚠ High flow mode should be activated to enjoy maximum Flow rate.

Once pressure set is reached, the insufflator stops insufflating. Insufflation restarts as soon as intracavitory pressure drops below pressure set.

⚠ Control the tightness of trocars where the gas input.

G. OVER-PRESSIONS

As soon as intracavitory pressure exceeds pressure set by 2 mmHg:

- The message " OVER PRESSURE" is displayed on the screen
- There is also a visual warning light [B]


After one second the insufflators open the external pressure relief valve to decrease the intracavitory pressure.

As soon as intracavitory pressure exceeds pressure set by 5 mmHg :

- The message " OVER PRESSURE" is displayed on the screen
 - There is also a visual warning light [B] and sound signal.
- After 3 seconds the device discharges the gas by opening the external pressure relief valve.


H. MONITORING THE FLOW

The instantaneous flow rate is indicated on the LCD screen in liters/minute.

I. MONITORING THE VOLUME OF GAS USED

The volume of gas used is shown in litres on the LCD display to the nearest deciliter


This meter is started by pressing on the button [F] and only resets to zero once the device is powered off.

- i** If necessary set it to zero via the menu.
- ⚠** Use only with medical grade CO₂ for coelioscopy complying with European standards.

J. MONITORING THE VOLUME OF GAS AVAILABLE

1. With gas cylinder

As soon as the cylinder pressure reaches 40 bars a updated cylinder icon is displayed, is flickering.
At 20 bars, a yellow light alarm light [B] and an alarm message "CYLINDER PRESSURE 20 bar" will be displayed.


- ⚠** From 10 bars, it is no longer possible to start the device (RUN/STOP key disabled); the cylinder must be changed.
- ⚠** Always have a new cylinder in the OR room.

2. With central gas wall CO₂ supply

As soon as the central gas wall CO₂ pressure reaches 2 bars, the insufflators inform the user via :

- A message on the display in the warning area "input pressure".
- A yellow light displayed.

[i] From, 2 bars insufflator is not allowed to re-start.

[i] From, 10 bars insufflators will automatically switch in "bottle" mode for safety reasons.

⚠ Use only with central GAS WALL CO₂ supply complying with European standards.

K. END OF OPERATION

Stop insufflation by pressing [F]; the STOP icon appears.

Disconnect the tubing from the trocar and the insufflator immediately to prevent any liquid or gas from flowing back to the device.

⚠ The used tubing must be disposed of in an appropriate container after use.

[i] The total volume of gas insufflated during the operation appears; reset the gas meter to zero if another operation is planned (see chapter VII - C).

Close the CO₂ gas cylinder.

FAULTS AND WARNING SIGNALS

WARNING MESSAGE	POSSIBLE CAUSE	ACTION TO BE TAKEN
END OF CYLINDER	Cylinder pressure at 40 bar	Prepare a replacement cylinder
	Cylinder pressure at 20 bar	Change the cylinder
	Gas cylinder valve closed	Open the valve
	Insufflator high-pressure hose connector disconnected	Turn off the device, check CO ₂ intake and restart the device
INPUT PRESSURE	Wall CO ₂ gas network failing	Check that the pressure delivered by the wall CO ₂ gas network is superior to 2 bars
OVER-PRESSURIZATION	Sporadic action on pneumoperitoneum	<p>None.</p> <ul style="list-style-type: none"> If the pressure is higher of more than 2mmHg to the pressure set, the device performs a light external desufflation of 500 ms. If the pressure is higher of more than 5 mmHg to the pressure set, the device performs a light external desufflation during 1 s.
	Tubing kinked	Stretch the tubing
	Trocar valve closed during the insufflation.	Open the trocar
TUBING NOT THERE	Tubing is missing	Plug the tubing
	The tubing is wrongly connected (desufflation actuator)	Connect the tubing correctly on the [C] and [D]
CALIBRATION ERROR	The tubing is connected	Disconnect the tubing from insufflator's outlet connector for 2 seconds only
REMOVE TUBING	AUTOMATIC TEST not run	Remove the tubing for 2 seconds, then replace
	Fault in tubing detection	Send device for servicing
OBSTRUCTED TUBING / OVERPRESSURE	Angled tubing	None. Immediate desufflation after each insufflation.
	Trocar valve closed during insufflation	Immediate desufflation after each insufflation. Open the trocar or stop the device.
CO ₂ NOT DETECTED	Gas cylinder valve closed	Open the valve
	Insufflator/bottle high-pressure hose disconnected	Connect the hose correctly to the insufflator and to the bottle
	CO ₂ gas wall network failing	Check that the pressure delivered by the wall CO ₂ gas network is superior to 2 bars
	Connection flexible hose insufflator/ wall network.	Connect carefully the hose to the CO ₂ gas wall network plug and to the insufflator.
OBLIGATORY MAINTENANCE	Number of operations performed before exceeded maintenance	Send device for Servicing
 SERVICE	Maintenance action essential	Send device for Servicing
NO POWER	Fuses damaged	Power switch off. Check and replace the 2.5A UR fuses if necessary. Switch it on.

 For any other problem, please contact your nearest after-sales service team

CLEANING / SERVICING / MAINTENANCE

A. CLEANING / DECONTAMINATION

⚠ Only instructions drawn up by the people responsible for sterilization in each hospital or health centre must be followed; these instructions will take precedence over the information contained in this manual, which is provided for guidance only.

⚠ Always disconnect the device's electrical socket before cleaning it.

After each use:

- discard disposable sterile tubing; do not try to sterilize it
- clean up any possible splashes of liquid on the insufflator by wiping them with a slightly dampened cloth.

⚠ The device must without fail be decontaminated before it is sent to after-sales service.

B. SERVICING / MAINTENANCE

⚠ Before 2000 surgeries or before one year of use a maintenance visit will need to be made by the manufacturer alone to ensure efficient and completely safe functioning.

This maintenance covers all the end-of-production check / inspection points and adjustments, checking of flow measuring system and pressure, inspection and replacement of the CO₂ internal filter.

SYMBOLS USED

DESCRIPTION	CORRESPONDING
Class 1 device	Class 1
Type BF device	
Warning	
UR time-lag fuse (conforms to UL standard)	
Record of equipotential	
Connector for CO ₂ port	
Date of manufacture MM/YYYY	
Manufacturer	
Complies with European directive 93/42/CEE	
Complies with the European directive 2002/96/CE on Waste Electrical and Electronic Equipment (DEEE or WEEE).	
Please read the manual !	

TECHNICAL SPECIFICATIONS

Pneumatics:

- CO₂ gas supply system: US 7/16" connector
- Pressure range :
 - CO₂ bottle : 10 to 60 Bars
 - Central gas wall supply : 3,5 to 4 Bars
- Maximum flow without loss of head (to the nearest decilitre):
 - 20 to 45 l/min in high flow (depending on the model)
 - 2 l/min in low flow
- Pressure set: 0 to 25 mmHg (Accuracy: 1mmHg)
- External desufflation valve
- Option of automatic low/high flow switching

Interface:

- 5.7-inch LCD screen showing: immediate flow, intra-cavity pressure, blow-off pressure, total volume of CO₂ used and low levels in the cylinder
- Thumbwheel navigation and validation system

Power supplies:

- Types of power supply: 100-230 V AC 50-60 Hz
- Protection by fuses: 2 x 2.5 AT - 250V time-lag fuses - UR
- Power used: 75VA

Mechanics:

- Dimensions: L=300mm, D=340mm, H=145mm
- Weight: 5 kg

Operating - transportation and Storage Environment:

- Operating temperature range: +10°C / +40°C
- Operating relative humidity range: 30 to 75 %
- Transport and storage temperature range: +10°C / +40°C
- Transport and storage relative humidity range: 20% to 85%
- Operating, transport and storage atmospheric pressure range: 700hPa to 1060hPa

Standards:

- Electrical protection: class 1 type BF (Note : Applied part is the end of the tubing)
- Conforms to standard IEC 60 601-1; with US and Canadian deviations
- Not protected against water (IPX0)
- Not adapted for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide

Delivered accessories:

- Main power cable
- Open-ended spanner intended for CO₂ high-pressure feed hose connector.
- Extended CO₂ filter for high pressure hose.
- User manual

ELECTROMAGNETIC COMPATIBILITY

12.1 Guide and declaration by manufacturer - electromagnetic emissions

This insufflator is designed to be used within the electromagnetic environment specified below. The user should ensure that it is indeed used within this environment

Emissions test	Compliance	Electromagnetic environment - Guide
RF emissions CISPR 11	Group 1	This insufflator only uses radioelectrical energy for its internal subsystems. Therefore, it emits very low RF energy and is not likely to interfere with nearby electronic devices.
RF emissions CISPR 11	Class A	
Harmonic emissions EN 61000-3-2	Conform	
Voltage fluctuations / flicker EN 61000-3-3	Conform	

12.2 Guide and declaration by the manufacturer - electromagnetic immunity

This insufflator is designed to be used within the electromagnetic environment specified below. The user should ensure that it is indeed used within this environment

Immunity test	CEI 60601 severity level	Compliance level	Electromagnetic environment - Guide
Electrostatic discharges EN 61000-4-2	± 6 kV via contact ± 8 kV via air	± 6 kV ± 8 kV	The floor should be in timber, concrete or tiling. If the floor is covered with a synthetic material, the relative humidity should be at least 30%.
Quick transients in bursts EN 61000-4-4	± 2 kV for power lines ± 1 kV for input/output lines	± 2 kV ± 1 kV	The quality of the main power supply should be that of a typical commercial or hospital environment.
Voltage shocks EN 61000-4-5	Differential mode ± 1 kV Common mode ± 2 kV	± 1 kV ± 2 kV	The quality of the main power supply should be that of a typical commercial or hospital environment.
Brown-outs, brief power failures and voltage variations EN 61000-4-11	• <5% U _r - for 10 ms • 40% U _r - for 100 ms • 70% U _r - for 500 ms • <5% U _r - for 5 s	<5% U _r 10 ms <40% U _r 100 ms <70% U _r 500 ms <5% U _r 5 s	The quality of the main power supply should be that of a typical commercial or hospital environment. If the user of this insufflator requires being able to continue to work during interruptions in the main power supply, it is recommended that this insufflator be powered via a UPS or battery.
Magnetic field at network frequency (50/60 Hz)	3 A/m	3 A/m	The magnetic field at the network frequency should be at a level characteristic of a location within a typical commercial or hospital environment.

Note: U_r is the nominal value of the power voltage applied during the test.

Guide and declaration by the manufacturer - electromagnetic immunity

This insufflator is designed to be used within the electromagnetic environment specified below. The user should ensure that it is indeed used within this environment

Immunity test	CEI 60601 severity level	Compliance level	Electromagnetic environment - Guide
Conducted RF EN 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	<p>Portable and mobile RF communication devices must not be used at a distance from this insufflator - including cables - that is less than the recommended separation distance calculated using applicable formulae according to the transmitter's frequency.</p> <p>Recommended separation distance</p> $d = 1,16 \sqrt{P}$ $d = 1,16 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,33 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>In which P is the maximum output power of the transmitter, in Watts (W), assigned by the manufacturer of the transmitter and the recommended separation distance (d) in meters (m).</p> <p>The field levels emitted by immobile RF transmitters - which are to be established by an electromagnetic measurement at the site - must be lower than the compliance level in each frequency band.</p> <p>Interference can occur nearby devices bearing the following symbol:</p> 
Radiated RF EN 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	

Note 1 At 80 MHz and 800 MHz, the higher frequency band applies.

Note 2 These recommendations may not apply in all situations. The propagation of electromagnetic waves is modified by absorption and reflection due to structures, objects and people.

a The field level of immobile transmitters, such as radio telephone base stations (cellular and cordless) and mobile terrestrial radio systems, amateur radio systems, AM/FM radio communication systems, and TV systems cannot be evaluated theoretically with precision. To ascertain the electromagnetic environment due to immobile RF transmitters, a site measurement must be performed. If a field level measured within the environment in which this insufflator is used exceeds the above applicable compliance levels, check that this insufflator operates in a satisfactory manner. If abnormal operation is observed, additional measurements should be taken, such as reorientation or relocation of the reference system.

b Outside the frequency band of 150 kHz to 80 MHz, the field level should be less than 3 V/m

12.3 Separation distances recommended between portable and mobile RF communications systems and this insufflator

This insufflator is designed to be used within an electromagnetic environment in which radiated RF interference is controlled. The user of this insufflator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication systems (transmitters) and this insufflator, such as recommended below, as a function of the maximum output power of the communications system.

Maximum assigned output power of the transmitter W	Separation distance as a function of the transmitter's frequency M		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1,16\sqrt{P}$	$d = 1,16\sqrt{P}$	$d = 2,33\sqrt{P}$
0.01	0.116	0.116	0.233
0.1	0.366	0.366	0.736
1	1.16	1.16	2.33
10	3.66	3.66	7.36
100	11.6	11.6	23.3

Note 1 At 80 MHz and 800 MHz, the higher frequency band applies.

Note 2 These recommendations may not apply in all situations. The propagation of electromagnetic waves is modified by absorption and reflection due to structures, objects and people.

For transmitters of which the maximum output power is not listed above, the recommended separation distance d in meters (m) can be established by using the equation applicable for the transmitter frequency, in which P is the transmitter's maximum output power in Watts (W) assigned by the transmitter's manufacturer.

Note 1 At 80 MHz and 800 MHz, the separation distance stated for the higher frequency range should be applied.

Note 2 These recommendations may not apply in all situations. The propagation of electromagnetic waves is modified by absorption and reflection due to structures, objects and people.

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