SOPRO
397

Full HD
1080P

USER MANUAL

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

SOPRO
COMEG

Sopro 397:

USER MANUAL	P.2
MANUEL D'UTILISATION	P.20
MANUALE PER UTILIZZATORE	P.38
MANUAL DE USUARIO	P.56
BETRIEBSANLEITUNG	P.74

CONTENTS

1	INTRODUCTION	3
2	PRESENTATION OF THE CAMERA	3
3	PRECAUTIONS FOR USE	4
4	REGULATORY RECOMMENDATION	5
	4.1 CONFORMITY	
	4.2 ELECTROMAGNETIC INTERFERENCE & ELECTROSTATIC DISCHARGES	
	4.3 MATERIO-VIGILANCE	
	4.4 END OF LIFE	
5	INSTALLATION	6
	5.1 CONNEXION OF THE BACK PANEL	
	5.2 50 / 60HZ	
	5.3 STARTING	
6	USER ADJUSTMENTS	7
	6.1 WHITE BALANCE OR AWB [S2]	
	6.2 FOCUS	
	6.3 FUNCTIONNING OF THE MENUS	
	6.4 APPLICATION FIELD MENU [S3]	
	6.5 VIDEO PARAMETERS MENU [S4]	
	6.6 OPTION MENU [S5]	
	6.7 TEXT GENERATOR [S9]	
7	DESCRIPTION OF THE BACK PANEL	10
	7.1 MAINS PLUG	
	7.2 VIDEO OUTPUTS	
	7.3 PS2-KEYBOARD	
	7.4 PERIPHERAL OUTPUTS	
	7.5 IDENTIFICATION AND SYMBOLS ON THE BACK PANEL	
8	RECOMMENDED CLEANING PROCEDURE	11
9	CLEANING THE OPTICAL SURFACE	12
10	AFTER SALES SERVICE	13
	10.1 THE STAND-BY GREEN LED (L1) DOES NOT LIGHT WHEN YOU POWER THE CAMERA	
	10.2 THE STAND-BY GREEN LED LIGHTS, BUT THERE IS NO IMAGE ON THE MONITOR	
	10.3 THE IMAGE IS FUZZY, UNIFORMLY WHITE	
	10.4 THE IMAGE IS TOO CLEAR OR TOO DARK	
11	TECHNICAL SPECIFICATIONS	15
	11.1 CAMERA CONTROL UNIT (CCU)	
	11.2 CAMERA HEAD	
12	ELECTROMAGNETIC COMPATIBILITY	17

INTRODUCTION

Thank you for your confidence in us proven by the purchase of this device.

To take full advantage and comply with all necessary precautions, you should carefully read this manual and be well acquainted with its contents.

Sentences containing the symbol  are points requiring special attention.

Sentences containing the symbol  are for information.

To facilitate the installation and the use of the camera we wanted to make his handbook more practical. Thus references in the page of presentation in the back of the cover (in the form of [D1] for example) are proposed for you to visualize easily the parts of the product concerned.

PRESENTATION OF THE CAMERA

This camera is a 3CCD Full High definition colour camera, 100% digital, with deported electronic. Its ergonomics, its automatic shutter, its good sensitivity and excellent resolution as well as its exceptional fidelity of colours make of it the ideal medical tool, especially for flexible and rigid endoscopies.

It includes:

- Camera Control Unit (CCU)
- Camera Head
- Power cable
- Video cables Peripheral cable
- User manual

 For the USA and Canada, you must use the “hospital grade” cable supplied and connect it to an “hospital grade” mains outlet.

This equipment is delivered in a cardboard box, to be kept for further transport.

Mandatory Accessories:

- 24mm HD objective (S_097_0001)
- or
- 16-34 optical Zoom (S_413024)

Recommended Accessories (as an option):

- PS2 keyboard
- A XENON light source 300 W
- A screen with a resolution of 1920 x 1080 pixels and a DVI input

PRECAUTION FOR USE

- ✓ Respect the use and storage conditions.
- ✓ Install the CCU in a clean, dry and well ventilated place.
- ✓ The installation of the equipment must be exclusively done by a competent technician trained by the manufacturer.
- ✓ To avoid risk of electrical shocks, fire, short-circuits or dangerous emanations, do not insert any metallic object in the apparatus.
- ✓ Do not expose the CCU to water and do not store it in damp places (risk of electric shock).
- ✓ Use only the accessories supplied with the apparatus or proposed as optional by the manufacturer.
- ✓ Do not use this camera near any flammable anesthetics.
- ✓ This camera must not be used in an ionizing environment.
- ✓ Do not expose this device to LASER rays
- ✓ This camera is not a sterile device.
- ✓ Before any use, you should look for any rough surfaces, sharp edges, protrusions that could cause a safety hazard.

REGULATORY RECOMMENDATION

4.1 CONFORMITY

This device has been designed and manufactured by a company with a certified quality system. It complies with the requirements of the European Directive 93/42/EEC related to medical devices. As a result, it complies in particular with the ad hoc electric safety (IEC) and Electromagnetic Compatibility (EMC) standards.

4.2 ELECTROMAGNETIC INTERFERENCE AND ELECTROSTATIC DISCHARGES

Although this product complies with the EMC standards, it is possible under very special circumstances that it may disturb other devices or even be disturbed in turn by other devices or by an unfavourable electromagnetic environment.

To prevent such situations, it is recommended:

- To take care of the quality of the electric network (and particularly the grounding system of all the devices and carriages) ;
- To keep the device away from ElectroMagnetic sources (e.g., compressor, motor, transformer, HF generator, etc.).

4.3 MATERIO VIGILANCE

Like for any medical system, this device is subject to the provisions pertaining to a materio vigilance; any serious dysfunction must therefore be promptly notified to the competent authorities and the manufacturer with the greatest possible precision.

Address of the manufacturer: see on the last page

Caution: Federal law restricts this device to sale by or on the order of a Physician

 This camera is intended to be used by qualified physicians in general and plastic surgery to provide access, illumination and allow observation or manipulation of body cavities, hollow organs, and canals.

4.4 END OF LIFE

This appliance is marked according to the European Directive 2002/96/EEC on Waste Electrical and Electronic Equipment (WEEE).

By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol  on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for recycling of electrical and electronic equipment.

Disposal must be carried out in accordance with local environmental regulations for waste disposal.

For more detailed information about treatment, recovery and recycling of this product, please contact your local city office your household waste disposal service or the dealer to whom you purchased the product from.

INSTALLATION

Install the control unit on your video furniture on a table and take the different accessories necessary.

5.1 CONNEXION OF THE BACK PANEL

- Connect the power supply cable to the mains plug of the camera [C2] and a video cable to one of the video outputs: select DVI [C3], it gives a better result than the others.
- Plug the other end of this video cable on the video input of your monitor.

5.2 50/60HZ

For SD outputs [C6, C7] you need to select according to the country the frequency switch position [S13]. Use this Switch when the camera is powered off.

5.3 STARTING

 To optimize the quality of the images, always connect the sensor prior to power the camera.

- Plug the connector of the sensor on the front panel of the control unit [C1].
- Check that a 24mm objective is correctly screwed on the C-Mount.
- Power the monitor.
- Power the camera thanks to switch located on the rear panel [S14].

The camera is on stand-by, press on the button [S1] to power the camera, the green led [L1] will blink quickly during few seconds then become stable, the picture appears.

Another press on [S1] will put the camera on Stand-by position.

 Before placing your telescope (or fiberscope) in front of the lens of the camera, you must check the condition of the telescope (clean distal lens, correct focus if necessary).

USER ADJUSTMENTS

6.1 WHITE BALANCE OR AWB [S2]

It can be performed by two means.

1. From the control unit, select the White Balance icon and validate (S2).
2. From one of buttons on the sensors according to your programming (Long press action).

You must proceed as follows:

- Once the camera is connected to the endoscope and the light source powered, film a white surface;
- Start the white balance, the message WHITE BALANCING appears on the screen.
- Keep filming a white surface until the message WHITE BALANCE OK is displayed.

6.2 FOCUS

Once the endoscope is connected and the light source powered, turn slowly the objective ring to find a position where objects are clear.

i To obtain a satisfying depth of field for the operation, make a distant focus; you will thus avoid making it too often.

6.3 FUNCTIONNING OF THE MENUS

A MENU is displayed directly on the screen. Press one of these buttons (S3) (S4) (S5) and a first Line appears. Select the parameter with the two buttons (S6) then validate with (S7) button or cancel with (S8).

i Parameters are automatically saved when validated

6.4 APPLICATION FIELD MENU [S3]

Select the specialty amongst the following sub-menu; it will automatically set the video parameters advised for the application:

LAPAROSCOPY 10mm
 LAPAROSCOPY 5mm
 ARTHROSCOPY/ENT 4mm
 ARTHROSCOPY/ENT 3mm
 HYSTERO/CYSTOSCOPY 4mm
 HYSTERO/CYSTOSCOPY 3mm
 FLEXIBLE
 User 1
 User2

6.5 VIDEO PARAMETERS MENU [S4]

The user can modify different parameters.

BRIGHTNESS	Change the luminosity of the picture
BLUE GAIN	Change colorimetry
RED GAIN	Changer colorimetry
COLORS	Change the saturation of the colors
Noise Reduction	Limit the noise in the picture
WINDOW	Shutter window analysis: surface of the image on which the camera analyses the quantity of light received, in order to adjust the speed of the shutter to obtain a right exposure time. To be changed according to the scope size
GAMMA	Change the luminosity between dark and bright areas
DETAILS	Change the sharpness of the picture to get 3D effect or at the opposite smooth the picture with the Antimoire filter for flexible scopes.
GAIN	Change the level of luminosity in low light condition

6.6 OPTION MENU [S5]

Languages	FR: French; IT: Italian; ES: Spanish; EN: English; DE: German																				
BUTTON TOP BUTTON LEFT BUTTON RIGHT	Following functions are available on each button: <table><tr><td>AWB</td><td>ZOOM -</td><td>P1/P2</td><td>DET. +</td><td>GAIN -</td></tr><tr><td>MIRROR</td><td>ZOOM</td><td>BRI. +</td><td>DET. -</td><td>GAIN</td></tr><tr><td>FREEZE</td><td>PER. 1</td><td>BRI. -</td><td>DETAILS</td><td>INFO</td></tr><tr><td>ZOOM +</td><td>PER. 2</td><td>BRIGHT</td><td>GAIN +</td><td></td></tr></table>	AWB	ZOOM -	P1/P2	DET. +	GAIN -	MIRROR	ZOOM	BRI. +	DET. -	GAIN	FREEZE	PER. 1	BRI. -	DETAILS	INFO	ZOOM +	PER. 2	BRIGHT	GAIN +	
AWB	ZOOM -	P1/P2	DET. +	GAIN -																	
MIRROR	ZOOM	BRI. +	DET. -	GAIN																	
FREEZE	PER. 1	BRI. -	DETAILS	INFO																	
ZOOM +	PER. 2	BRIGHT	GAIN +																		
Factory Settings	Re-setting of the camera with factory parameters																				

6.7 TEXT GENERATOR [S9]

- Connect a keyboard on the rear panel [C8]
- Press on [S9] to make the text appear
- Write the text then it appears on one line

i If the text is not displayed on the screen it does not appear on the video signal output.

DESCRIPTION OF THE BACK PANEL

7.1 MAINS PLUG

The electric power supply of the camera is done through a plug [C2] which must be linked to the local supply circuit thanks to the electric cable delivered with the system. This plug includes a fuse-trap [F1] and the general power switch. The fuses used are of type T2.5A - 250V.

During a changing of fuse, it is peremptory to disconnect the camera from the mains and to use a fuse of the same kind. The T of “T2.5A” means “temporised”. Use only fuses marked UR.

7.2 VIDEO OUTPUTS

The following video outputs are available on this camera: DVI [C3], HD-SDI [C4], RGBS [C5], Y/C [C6], and composite [C7].

 DVI output is the connection which ensures the best picture quality with 1080p format on short distances.

 HD-SDI outputs allow carrying the video signal on long distances (deported screen up to 100m....)

7.3 PS2-KEYBOARD

Connect the PS2 certified keyboard (delivered as an option) on the dedicated connector [C8].

7.4 PERIPHERAL OUTPUTS

Two peripheral outputs [C9] are available and allow through a jack cable 3,5mm to drive 2 peripherals (ex : a printer or video recorder) from the buttons located on the sensor [S10,11,12].

 These outputs must be configured though the OPTION menu.

 The devices which are connected to the “video out” and “peripheral” plugs must comply with the CEI 950 standard.

7.5 IDENTIFICATION AND SYMBOLS ON THE BACK PANEL

The indications and symbols mentioned on the back panel allow identifying the camera in conformity with the International standards IEC 60601-1, IEC 60601-2-18, IEC 60417 and ISO 980 (See Annex 1).

RECOMMENDED CLEANING PROCEDURE

⚠ Whatever the procedure followed for the decontamination is, the use of a sterile sleeve is compulsory during all the duration of the surgical operation.

The sensor and the objective of the camera are designed for a cold disinfection by immersion of glutaraldehyde type (Lab. Anios). All other methods of disinfection are prohibited. The damages due to these other methods could not be taken over by the manufacturer.

The camera connector equipped with its waterproofness cap can be immersed. It is then peremptory to dry carefully the interior of the cap before re-use. The traces of water that must be dried appear when you remove the cap from the sensor (trickling and condensation).

We advise you to immerse the sensor with its objective screwed on it, in order to avoid having stains on the glass before the CCD.

It is peremptory to respect the pre disinfection procedure before the disinfection one and, after disinfection, to rinse carefully with sterile water all the parts which have been in contact with the disinfectant.

Cleaning - pre disinfection	Disinfection	Sterilization
Immersion into a cleaning solution of enzymatic type (Aniosyme DD1, Hexanios G+R) or cleaning of the surface (Anios sweepers).	Immersion into a glutaraldehyde solution (2% or 3.2%) (Steranios 2%)	Sterilization by process STERRAD® 100S

⚠ You must select one of the two solutions indicated for this type of equipment and use it exclusively during all the lifetime of the product for the pre-disinfection and disinfection processes.

⚠ This camera is not autoclavable

⚠ This camera is not compatible with automatic washers-disinfectors.

⚠ Existing alkaline solution for the pre-disinfection of some medical devices are FORBIDDEN for the pre-disinfection of our cameras.

⚠ It is peremptory to thoroughly rinse with sterile water the parts which have been in contact with the disinfectant.

⚠ Use unweaved compresses for the drying of the distal lens to avoid scratching it.

⚠ The procedures described in this chapter are given to you as advice. They can in no case substitute for official recommendations or directives.

⚠ Decontamination being linked to the products, methods and/or selected tools, it remains under the entire responsibility of the people in charge.

CLEANING THE OPTICAL SURFACE

Observation can be impeded by a dirty lens surface. There should not be any fingerprints or marks on the surface of the lens. To avoid scratching the lens surface, never use an abrasive cloth or sponge to clean it.

To remove dust, dirt and other impurities which do not come from the patient, wipe using a soft, lint-free cloth soaked in ethanol or isopropyl alcohol at 70%. Wear gloves which protect against chemical products during cleaning.

Always ensure that the equipment is completely dry before using it.

 Use unweaved compresses for the drying of the distal lens to avoid scratching it.

AFTER SALE SERVICE

No special maintenance is necessary on the camera. However, here is a list of common problems which are easy to troubleshoot. For all other cases, contact your Maintenance Service.

 We remind you that the mishandling is not guaranteed.

10.1 THE STAND-BY GREEN LED (L1) DOES NOT LIGHT WHEN YOU POWER THE CAMERA

- Check if the mains plug on the rear panel of the camera is well connected to the mains and the general switch is ON ("I").
- Check the good state of the fuses (use only fuse of the kind of those specified on the rear panel).

10.2 THE STAND-BY GREEN LED LIGHTS, BUT THERE IS NO IMAGE ON THE MONITOR

- Check if the sensor is well connected to the front panel (in the opposite case, a colour bar will be displayed on the monitor) and that it films a workable image.
- Check if the camera is well connected to the monitor (good cable, correctly plugged).
- If your video chain holds several elements, try to connect directly the camera to the monitor in order to validate that the problem is not located at the level of one of your peripherals.
- Check if the monitor is well powered, switched on the good video input and if the images adjustments are in a middle position (colour, light and contrast).
- Check your light source, your light cable and your telescope.

10.3 THE IMAGE IS FUZZY, UNIFORMLY WHITE

- Check if there is no blur on the lens.
- Check if the sensor does not aim at a too lighted object.
- Check the focus of the lens.

10.4 THE IMAGE IS TOO CLEAR OR TOO DARK.

- Check if the “BRIGHTNESS” parameter is not at its maximum or minimum.
- Check if the “WINDOW” parameter used is not too small or too big.
- Check if the camera head is correctly connected to the control unit.
- Turn off and power again the camera (no Hot plug of the camera head).

If the problem remains and if you are led to send back the camera, once disinfected, be careful to do it in its original wrapping.

Likewise, you must send us the whole system (control unit, sensor, lenses, and cables). Be kind enough to join an explanatory note relative to the observed default to the expedition bill.

 **The device must be disinfected before being sent back to the manufacturer.**

At the reception of your equipment, check its state and enter a reservation on the delivery bill if necessary, then confirm them by recommended letter addressed to the carrier within 48 hours. In the case where equipment shipped by us suffers damages during the travel, the amount of the repairing will be charged whether to the carrier if the reservation were entered in the delays, or to the consignee in the contrary.

TECHNICAL SPECIFICATIONS

11.1 CAMERA CONTROL UNIT (CCU)

Main Functions

- Class I product
- 3 Programmable buttons
- DIGITAL ZOOM Up to 2.5x
- FREEZE & MIRROR
- MENU by Application
- Digital processing 3 x 19 bits
- Automatic white balance, Range: 2800K to 6500K
- Two peripherals triggering
- 5 languages (English, French, German, Spanish, Italian)
- Color bar generator
- Automatic and manual adjustment of color
- Color Matrix 12 axes
- Automatic saving of the adjustments parameters
- Electronic shutter : OFF / Manual / ELC (AREA setting, PEAK/AVERAGE setting)
- Automatic (1/50 to 1/100 000)
- Adjustable shutter area
- Adjustments displayed on the monitor
- Text generator with Keyboard connexion
- Gamma adjustment
- Automatic and manual adjustment of light intensity
- Automatic adjustment of sensitivity to low light (auto exposure control)
- Automatic black level adjustment
- Details adjustment
- Digital Noise Reduction: OFF / LOW / HIGH
- Automatic compensation function for the blemish pixels

Video outputs

- 1 VB5: 1.0 V [P-P] / 75 ohms Sync Negative,
- 1 Y/C: Y: 1.0 V [P-P] / 75 ohms, C: 0.286 V [P-P] (NTSC), C: 300 Vp-p V [P-P] (PAL) / 75 ohms Sync Negative,
- 1 RGB: 3x0.7 V [P-P] & 0.3V [P-P] / 75 ohms
- 2 HD-SDI: 0.8 V [P-P] / 75 ohms SMPTE292M
- 1 DVI-D 3x0,7Vcc/& 0,3Vcc/75Ω)

Environment

- Operating Temperature range: +10°C / +40°C
- Operating Relative Humidity range: 30 to 75 %
- Transport and storage Temperature range: -10°C / +45°C
- Transport and storage Relative Humidity range: 20 to 85 %
- Operating, Transport and Storage Atmospheric pressure range: 700hPa to 1060hPa
- Not protected against waterfalls (IPX0)
- Not adapted to a use close to a mixture of inflammable anaesthetic with air or oxygen or nitrous oxide.
- Comply with the European standard 93/42/CEE.
- Comply with the international standards IEC 60601-1; IEC 60601-2-18 ; IEC 60417 and EN 980
- Size of the control unit : W : 300mm; D : 311mm; H : 125 mm
- Weight of the control unit : 2500 g

Typical Operating Requirements:

- Input Voltage: 100 - 230 VAC @50/60 Hz
- Power Consumption: 45 VA
- 2 T 2.5A - 250 V fuses (Only fuses marked UR)
- Continuous Service

11.2 CAMERA HEAD

- BF Type applied part
- Sensor : 3 CCD 1/3"
- Capturing aspect ratio on the CCD: 16 / 9 ("wide - screen format acquisition")
- Progressive Scan Technology
- 1080p resolution HDTV format: Resolution : 1920*1080 pixels
- 14 bits processing
- Definition : 1080 lines
- Sensitivity : 1,5 lux
- Signal / Noise ratio : 70dB
- Cable length : 3 meters
- C-Mount
- Dimensions of the sensor: 37*50*95 mm
- Weight of the sensor: 250 g

ELECTROMAGNETIC COMPATIBILITY

12.1 GUIDE AND DECLARATION BY MANUFACTURER - ELECTROMAGNETIC EMISSIONS

The reference system, "3 ccd camera" is designed to be used within the electromagnetic environment specified below. The user should ensure that it is indeed used within this environment

Emissions test	Compliance	Electromagnetic environment - Guide
RF emissions CISPR 11	Group 1	The reference system, " 3 ccd camera " only uses radioelectrical energy for its internal subsystems. Therefore, it emits very low RF energy and is not likely to interfere with nearby electronic devices.
RF emissions CISPR 11	Class A	The reference system, " 3 ccd camera ", must be used in all premises other than residential premises and premises connected directly to the public low-voltage power distribution network used for powering residential buildings.
Harmonic emissions EN 61000-3-2	Class A	
Voltage fluctuations / flicker EN 61000-3-3	Conform	

12.2 GUIDE AND DECLARATION BY THE MANUFACTURER - ELECTROMAGNETIC IMMUNITY

The reference system, " 3 ccd camera " is designed to be used within the electromagnetic environment specified below. The user should ensure that it is indeed used within this environment

Immunity test	CEI 60601 severity level	Compliance level	Electromagnetic environment - Guide
Electrostatic discharges EN 61000-4-2	± 6 kV via contact ± 8 kV via air	± 6 kV ± 8 kV	The floor should be in timber, concrete or tiling. If the floor is covered with a synthetic material, the relative humidity should be at least 30%.
Quick transients in bursts EN 61000-4-4	± 2 kV for power lines ± 1 kV for input/output lines	± 2 kV ± 1 kV	The quality of the main power supply should be that of a typical commercial or hospital environment.
Voltage shocks EN 61000-4-5	Differential mode ± 1 kV Common mode ± 2 kV	± 1 kV ± 2 kV	The quality of the main power supply should be that of a typical commercial or hospital environment.
Brown-outs, brief power failures and voltage variations EN 61000-4-11	<5% U _r - for 10 ms 40% U _r - for 100 ms 70% U _r - for 500 ms <5% U _r - for 5 s	<5% U _r 10 ms <40% U _r 100 ms <70% U _r 500 ms <5% U _r 5 s	The quality of the main power supply should be that of a typical commercial or hospital environment. If the user of the reference system, " 3 ccd camera ", requires to be able to continue to work during interruptions in the main power supply, it is recommended that the reference system, " 3 ccd camera " be powered via a UPS or battery.
Magnetic field at network frequency (50/60 Hz)	3 A/m	3 A/m	The magnetic field at the network frequency should be at a level characteristic of a location within a typical commercial or hospital environment.
Conducted RF EN 61000-4-6 Radiated RF EN 61000-4-3	3 V _{rms} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3V 3V/m	Portable and mobile RF communication devices must not be used at a distance from the reference system, " 3 ccd camera " - including cables - that is less than the recommended separation distance calculated using applicable formulae according to the transmitter's frequency. Recommended separation distance $d = 1,16\sqrt{P}$ $d = 1,16\sqrt{P}$ 80 MHz à 800 MHz $d = 2,33\sqrt{P}$ 800 MHz à 2,5 GHz In which P is the maximum output power of the transmitter, in Watts (W), assigned by the manufacturer of the transmitter and the recommended separation distance (d) in meters (m). The field levels emitted by immobile RF transmitters - which are to be established by an electromagnetic measurement at the site - must be lower than the compliance level in each frequency band. Interference can occur nearby devices bearing the following symbol: 

NOTE : U_T is the nominal value of the power voltage applied during the test.

NOTE 1 : At 80 MHz and 800 MHz, the higher frequency band applies.

NOTE 2 : These recommendations may not apply in all situations. The propagation of electromagnetic waves is modified by absorption and reflection due to structures, objects and people.

a The field level of immobile transmitters, such as radio telephone base stations (cellular and cordless) and mobile terrestrial radio systems, amateur radio systems, AM/FM radio communication systems, and TV systems cannot be evaluated theoretically with precision. To ascertain the electromagnetic environment due to immobile RF transmitters, a site measurement must be performed. If a field level measured within the environment in which the reference system, " 3 ccd camera ", is used exceeds the above applicable compliance levels, check that the reference system, " 3 ccd camera ", operates in a satisfactory manner. If abnormal operation is observed, additional measurements should be taken, such as reorientation or relocation of the reference system.

b Outside the frequency band of 150 kHz to 80 MHz, the field level should be less than 3 V/m

12.3 SEPARATION DISTANCES RECOMMENDED BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS SYSTEMS AND THE REFERENCE SYSTEM, "3 CCD CAMERA"

The reference system, "3 ccd camera", is designed to be used within an electromagnetic environment in which radiated RF interference is controlled. The user of the reference system, "3 ccd camera", can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication systems (transmitters) and the reference system, "3 ccd camera", such as recommended below, as a function of the maximum output power of the communications system.

Maximum assigned output power of the transmitter W	Separation distance as a function of the transmitter's frequency M		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1,16\sqrt{P}$	$d = 1,16\sqrt{P}$	$d = 2,33\sqrt{P}$
0.01	0.116	0.116	0.233
0.1	0.366	0.366	0.736
1	1.16	1.16	2.33
10	3.66	3.66	7.36
100	11.6	11.6	23.3

For transmitters of which the maximum output power is not listed above, the recommended separation distance d in meters (m) can be established by using the equation applicable for the transmitter frequency, in which P is the transmitter's maximum output power in Watts (W) assigned by the transmitter's manufacturer.

NOTE 1 : At 80 MHz and 800 MHz, the separation distance stated for the higher frequency range should be applied.

NOTE 2 : These recommendations may not apply in all situations. The propagation of electromagnetic waves is modified by absorption and reflection due to structures, objects and people.

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